


Axygen® Pipet Tips

High quality pipet tips for all
your liquid handling needs



AXYGEN®

A Corning Brand



Axygen® Pipet Tips

Overview	2
Axygen MultiRack Tip System.....	3
Universal Filter Tips	3
Universal Tips.....	3
Axygen LTS Compatible Tips	5
Axygen HybridRack™ Pipet Tips	6
Axygen Traditional Pipet Tips	7
Universal Filter Tips	7
Universal Tips.....	8
Related Products.....	8

Product Ordering Information

Products may not be available in all markets. For information on Purchasing Options, Terms and Conditions of Sale, Return and Repair Policies, and Warranty/Guarantee Registration, visit www.corning.com/how-to-buy.


Products may not be available in all markets.

Overview

Axygen® universal pipet tips are designed to fit a wide variety of single- and multi-channel pipettors. Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject.

The Axygen product offering now also includes the Axygen MultiRack System, which is our updated reload and racking system. In addition to the universal fit tips packaged in the MultiRack System, we now also offer LTS Compatible tips in the MultiRack System, and if you are worried about your carbon footprint, the Axygen HybridRack™ is a perfect option for you.

- ▶ **Automated assembly** – Ensures product is free from contamination.
- ▶ **ISO 13485 and ISO 14001** – Certified quality management systems supply stringent quality processes to ensure the highest quality finished good.
- ▶ **In-house testing** – In-house QA lab allows product to be validated for fit, form, and function prior to release.
- ▶ **Brand you can trust** – A long track record of innovation and value as one of the Corning trusted brands.
- ▶ **Broad product offering** – Corning offers many different styles of pipet tips, including racked, filtered, reloads, stack racks, low retention, wide bore, and more. With so many options, you are sure to find a tip that works for you.



Packaging Solutions




Bulk



Reload



Racked



Tip Station



HybridRack



Select products are included in Corning EcoChoice™. All products included in the program are produced, packaged, or distributed in an environmentally friendly manner following United States FTC Guidelines. To learn more, go to www.corning.com/EcoChoice.

Axygen® MultiRack Tip System



Engineered to fit most pipettors on the market, the Axxygen MultiRack Tip System is a complete pipet tip system offering in various packaging styles to meet every lab's needs. The MultiRack Tip System offers standard and extended-length sizes in volumes from 0.1 µL to 1250 µL. Racked and reload packaging styles feature preloaded tips in easy-load, color-coded deck inserts. With a small benchtop footprint and reusable rack, this system minimizes environmental impact and maximizes lab space. Axxygen MultiRack Tip System products are made in the USA and meet ISO standards.

Axygen MultiRack Universal Filter Pipet Tips

- ▶ Racked
- ▶ Sterile
- ▶ RNase-/DNase-free
- ▶ Human DNA-free
- ▶ PCR inhibitor-free
- ▶ Non-pyrogenic
- ▶ USP Class VI polypropylene



Standard Surface Cat. No.	Maxymum Recovery [®] Surface Cat. No.	Volume	Description	Quantity
MRF-10XT-R-S	MRF-10XT-L-R-S	10 µL	Extended length, graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
MRF-20-R-S	MRF-20-L-R-S	20 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
MRF-100-R-S	MRF-100-L-R-S	100 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
MRF-200NX-R-S	MRF-200NX-L-R-S	200 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
MRF-200XT-R-S	MRF-200XT-L-R-S	200 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case
MRF-300XT-R-S	MRF-300XT-L-R-S	300 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case
MRF-1000XT-R-S	MRF-1000XT-L-R-S	1250 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case

Axygen MultiRack Universal Pipet Tips


10 µL Extended Length, Graduated, Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Packaging	Sterile	Quantity
MR-10XT	-	Bulk	No	1,000/pack, 10 packs/case
MR-10XT-RL	-	Reload	No	96/deck, 10 decks/pack, 4 packs/case
MR-10XT-R	-	Racked	No	96/rack, 10 racks/pack, 5 packs/case
MR-10XT-R-S	MR-10XT-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
MR-10XT-TS	-	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
MR-10XT-TS-S	-	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

Axygen® MultiRack Universal Pipet Tips (continued)

200 µL, Graduated, Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Packaging	Sterile	Quantity
MR-200	-	Bulk	No	1,000/pack, 10 packs/case
MR-200-RL	-	Reload	No	96/deck, 10 decks/pack, 4 packs/case
MR-200-R	-	Racked	No	96/rack, 10 racks/pack, 5 packs/case
MR-200-R-S	MR-200-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
MR-200-TS	-	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
MR-200-TS-S	-	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

300 µL, Graduated, Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Packaging	Sterile	Quantity
MR-300-RL	-	Reload	No	96/deck, 10 decks/pack, 4 packs/case
MR-300-R	-	Racked	No	96/rack, 10 racks/pack, 5 packs/case
MR-300-R-S	MR-300-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

300 µL Graduated, Extended Length, Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Packaging	Sterile	Quantity
MR-300XT-R	-	Racked	No	96/rack, 8 racks/pack, 5 packs/case
MR-300XT-R-S	MR-300XT-L-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case

1250 µL Graduated, Extended Length, Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Packaging	Sterile	Quantity
MR-1000XT	-	Bulk	No	768/pack, 10 packs/case
MR-1000XT-RL	-	Reload	No	96/deck, 8 decks/pack, 4 packs/case
MR-1000XT-R	-	Racked	No	96/rack, 8 racks/pack, 5 packs/case
MR-1000XT-R-S	MR-1000XT-L-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case

Axygen® LTS Compatible Tips



Axygen LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch pipettors, all with the quality you have come to expect from Corning. By reducing the force required to seal and eject, Axxygen LTS Compatible tips can help you avoid using excess force to obtain an effective seal between the pipettor and pipet tip. Axxygen LTS Compatible tips are available in your favorite packaging formats including Axxygen HybridRack™, a sustainable packaging choice. All LTS tips are composed of USP Class VI polypropylene and manufactured in an ISO 13485 and ISO 14001 United States based facility.

- RNase-/DNase-free
- Human gDNA-free
- PCR inhibition-free
- Nonpyrogenic
- USP Class VI polypropylene



20 µL LTS Compatible Extended Length, Graduated, Fine Point

Cat. No.	Packaging	Sterile	Quantity
LTF-20-R-S	Racked, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
LTFER-20-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
LT-20-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
LT-20-R	Racked	No	96/rack, 10 racks/pack, 5 packs/case
LT-20-RL	Reload stack	No	96/deck, 10 decks/pack, 4 packs/case
LT-20-TS	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
LT-20-TS-S	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case



200 µL LTS Compatible, Graduated, Fine Point

Cat. No.	Packaging	Sterile	Quantity
LTF-200-R-S	Racked, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
LTFER-200-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
LT-200-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
LT-200-R	Racked	No	96/rack, 10 racks/pack, 5 packs/case
LT-200-RL	Reload stack	No	96/deck, 10 decks/pack, 4 packs/case
LT-200-TS	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
LT-200-TS-S	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case



1000 µL LTS Compatible, Graduated, Fine Point

Cat. No.	Packaging	Sterile	Quantity
LTF-1000-R-S	Racked, filtered	Yes	96/rack, 8 racks/pack, 5 packs/case
LTFER-1000-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
LT-1000-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case
LT-1000-R	Racked	No	96/rack, 8 racks/pack, 5 packs/case
LT-1000-RL	Reload stack	No	96/deck, 8 decks/pack, 4 packs/case

Axygen® HybridRack™ Pipet Tips



Born from a desire to reduce the amount of waste happening in laboratories, Axxygen HybridRack pipet tips are an example of our innovative spirit leading to a more sustainable future. By replacing the traditional plastic rack and lid normally found in racked pipet tip products with a paperboard rack and lid, we can reduce the amount of plastic per case by 70%, or approximately two pounds per case of 200 μ L tips. The HybridRack is strong enough to support multi-channel pipetting, and as soon as the rack is empty, it can be recycled using any curbside recycling program.

- ▶ Sterile
- ▶ RNase-/DNase-free
- ▶ Human DNA-free
- ▶ PCR inhibitor-free
- ▶ Non-pyrogenic
- ▶ USP Class VI polypropylene

10 μ L Extended Length, Graduated, Fine Point

Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Filtered	Sterile	Quantity
ERF-10XT-R-S	ERF-10XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
ER-10XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

20 μ L, 100 μ L, and 200 μ L Graduated, Fine Point

Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Filtered	Sterile	Quantity
ERF-20-R-S	ERF-20-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
ERF-100-R-S	ERF-100-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
ER-200-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

200 μ L and 300 μ L Graduated, Fine Point

Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Filtered	Sterile	Quantity
ERF-200NX-R-S	ERF-200NX-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
ER-300NX-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

300 μ L Graduated, Extended Length, Fine Point

Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Filtered	Sterile	Quantity
ERF-300XT-R-S	ERF-300XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
ER-300XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case


1250 μ L Graduated, Extended Length, Fine Point

Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Filtered	Sterile	Quantity
ERF-1000XT-R-S	ERF-1000XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
ER-1000XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

To learn more about the Corning HybridRack tips, see the Axxygen HybridRack Pipet Tips Infographic (CLS-AN-782) at www.corning.com/lifesciences.



Axygen® Traditional Pipet Tips



Known around the world for industry leading quality, accuracy, and precision, Axxygen Traditional pipet tips have been the preferred liquid handling solution for laboratories around the world for the last 20+ years. Axxygen Traditional pipet tips are designed to fit a wide variety of single- and multi-channel pipettors. Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject. Featuring tips with extended length, wide bore designs for viscous solutions, and our trademarked Maxymum Recovery surface, Axxygen Traditional pipet tips are sure to meet and exceed your liquid handling expectations.

- ▶ RNase-/DNase-free
- ▶ Human gDNA-free
- ▶ Nonpyrogenic
- ▶ USP Class VI polypropylene

Axygen Traditional Universal Filter Pipet Tips

- ▶ Racked
- ▶ Filtered
- ▶ Sterile



Standard Surface Cat. No.	Maxymum Recovery® Surface Cat. No.	Volume (µL)	Description	Quantity
TF-300-R-S	TF-300-L-R-S	0.5 - 10	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
TF-400-R-S	TF-400-L-R-S	0.5 - 10	Extended length, graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
TF-420-R-S	TF-420-L-R-S	0.5 - 20	Extended length, graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
TXLF-10-R-S	TXLF-10-L-R-S	0.5 - 10	Extended length, fine point	96/rack, 10 racks/pack, 5 packs/case
TF-20-R-S	TF-20-L-R-S	0.5 - 20	Beveled	96/rack, 10 racks/pack, 5 packs/case
TF-50-R-S	TF-50-L-R-S	0.5 - 50	Graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
TF-100-R-S	TF-100-L-R-S	1 - 100	Graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
TF-200-R-S	TF-200-L-R-S	1 - 200	Graduated, beveled	96/rack, 10 racks/pack, 5 packs/case
TF-205-WB-R-S	TF-205-WB-L-R-S	1 - 200	Wide bore	96/rack, 10 racks/pack, 5 packs/case
TF-350-R-S	TF-350-L-R-S	1 - 300	Fine point	96/rack, 10 racks/pack, 5 packs/case
TF-1000-R-S	TF-1000-L-R-S	100 - 1000	Beveled	100/rack, 10 racks/pack, 5 packs/case
TF-1005-WB-R-S	TF-1005-WB-L-R-S	1 - 1000	Wide bore	100/rack, 10 racks/pack, 5 packs/case

Axygen® Traditional Universal Pipet Tips

10 µL Graduated, Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-300	T-300-L	0.5 - 10	Bulk	No	1,000/pack, 20 packs/case
T-300-R	T-300-L-R	0.5 - 10	Racked	No	96/rack, 10 racks/pack, 5 packs/case
T-300-R-S	T-300-L-R-S	0.5 - 10	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
T-300-STK	T-300-L-STK	0.5 - 10	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case
T-300-STK-S	T-300-L-STK-S	0.5 - 10	Stack rack	Yes	96/layer, 5 layers/stack, 10 stacks/case
RFL-300-C	-	0.5 - 10	Reload	No	96 tips/deck, 10 decks/pack, 5 packs/case

10 µL Extended Length, Graduated, Beveled



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-400	T-400-L	0.5 - 20	Bulk	No	1,000/pack, 20 packs/case
T-400-R	T-400-L-R	0.5 - 20	Racked	No	96/rack, 10 racks/pack, 5 packs/case
T-400-R-S	T-400-L-R-S	0.5 - 20	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

10 µL Extended Length, Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
TXL-10	TXL-10-L	0.5 - 20	Bulk	No	1,000/pack, 20 packs/case
TXL-10-R	TXL-10-L-R	0.5 - 20	Racked	No	96/rack, 10 racks/pack, 5 packs/case
TXL-10-R-S	TXL-10-L-R-S	0.5 - 20	Racked	Sterile	96/rack, 10 racks/pack, 5 packs/case

200 µL Beveled



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-200-C	T-200-C-L	1 - 200	Bulk	No	1,000/pack, 20 packs/case
T-200-C-R	T-200-C-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
T-200-C-R-S	T-200-C-L-R-S	1 - 200	Racked	Sterile	96/rack, 10 racks/pack, 5 packs/case
T-200-C-STK	T-200-C-L-STK	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case

200 µL Beveled, Yellow



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-200-Y	-	1 - 200	Bulk	No	1,000/pack, 20 packs/case
T-200-Y-R	-	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
T-200-Y-R-S	-	1 - 200	Racked	Sterile	96/rack, 10 racks/pack, 5 packs/case
T-200-Y-STK	-	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case

Axygen® Traditional Universal Pipet Tips (continued)

200 µL Reference-compatible, Beveled, Graduated



Standard Surface Cat. No.	Maxymum Recovery® Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
TR-222-C	TR-222-C-L	1 - 200	Bulk	No	1,000/pack, 20 packs/case
TR-222-C-R	TR-222-C-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
TR-222-C-R-S	TR-222-C-L-R-S	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
TR-222-C-STK	TR-222-C-L-STK	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case
TR-222-C-STK-S	TR-222-C-L-STK-S	1 - 200	Stack rack	Yes	96/layer, 5 layers/stack, 10 stacks/case

200 µL Reference-compatible, Beveled, Graduated, Yellow



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
TR-222-Y	-	1 - 200	Bulk	No	1000/pack, 20 packs/case
TR-222-Y-R	-	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
TR-222-Y-R-S	-	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
TR-222-Y-STK	-	1 - 200	Stack rack	No	96/layer, 5 layers/stack, 10 stacks/case
RFL-222-C	-	1 - 200	Reload	No	96 tips/deck, 10 decks/pack, 5 packs/case
RFL-R	-	-	Empty Rack for RFL-300-C and RFL-222-C	No	10 racks/pack, 5 packs/case

200 µL Wide Bore



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-205-WB-C	-	1 - 200	Bulk	No	1,000/pack, 20 packs/case
T-205-WB-C-R	-	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
T-205-WB-C-R-S	-	1 - 200	Racked	Sterile	96/rack, 10 racks/pack, 5 packs/case

300 µL Fine Point



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-350-C	T-350-C-L	1 - 300	Bulk	No	1,000/pack, 10 packs/case
T-350-C-R	T-350-C-L-R	1 - 300	Racked	No	96/rack, 10 racks/pack, 5 packs/case
T-350-C-R-S	T-350-C-L-R-S	1 - 300	Racked	Sterile	96/rack, 10 racks/pack, 5 packs/case

1000 µL Beveled, Clear



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-1000-C	T-1000-C-L	100 - 1000	Bulk	No	1,000/pack, 5 packs/case
T-1000-C-R	T-1000-C-L-R	100 - 1000	Racked	No	100/rack, 10 racks/pack, 5 packs/case
T-1000-C-R-S	T-1000-C-L-R-S	100 - 1000	Racked	Yes	100/rack, 10 racks/pack, 5 packs/case

Axygen® Traditional Universal Pipet Tips (continued)

1000 µL Beveled, Blue



Standard Surface Cat. No.	Maxymum Recovery® Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-1000-B	-	100 - 1000	Bulk	No	1,000/pack, 5 packs/case
T-1000-B-R	-	100 - 1000	Racked	No	100/rack, 10 racks/pack, 5 packs/case
T-1000-B-R-S	-	100 - 1000	Racked	Yes	100/rack, 10 racks/pack, 5 packs/case
RFL-1200-C	-	1 - 1200	Reload	No	96 tips/deck, 8 decks/pack, 6 packs/case
RFL-R-1200	-	-	Empty Rack for RFL-1200-C	No	10 racks/pack, 3 packs/case


1000 µL Wide Bore



Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Packaging	Sterile	Quantity
T-1005-WB-C	-	100 - 1000	Bulk	No	1,000/pack, 5 packs/case
T-1005-WB-C-R	T-1005-WB-C-L-R	100 - 1000	Racked	No	100/rack, 10 racks/pack, 5 packs/case
T-1005-WB-C-R-S	-	100 - 1000	Racked	Yes	100/rack, 10 racks/pack, 5 packs/case

Related Products

Axygen SnapLock Microcentrifuge Tubes



Cat. No.	Description	Quantity
MCT-060-C	0.6 mL Snaplock microcentrifuge tubes, MaxyClear	1,000/pack, 10 packs/case
MCT-060-C-S	0.6 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	100/bag, 5 bags/pack, 10 packs/case
MCT-150-C	1.5 mL Snaplock microcentrifuge tubes, MaxyClear	5,000/pack, 10 packs/case
MCT-150-C-S	1.5 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	50/bag, 5 bags/pack, 10 packs/case
MCT-175-C	1.7 mL Snaplock microcentrifuge tubes, MaxyClear	500/pack, 10 packs/case
MCT-175-C-S	1.7 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	250/pack, 10 packs/case
MCT-200-C	2.0 mL Snaplock microcentrifuge tubes, MaxyClear	5,000/pack, 10 packs/case
MCT-200-C-S	2.0 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	50/bag, 5 bags/pack, 10 packs/case
MCT-500-C	5.0 mL Snaplock microcentrifuge tubes, MaxyClear	250/pack, 5 packs/case
MCT-500-C-S	5.0 mL Snaplock microcentrifuge tubes, MaxyClear, sterile	250/pack, 5 packs/case

NOTE: Other options are available, including full cases of different colors, assorted colors in the same case, capless tubes, and tubes with Axygen Maxymum Recovery surface.

For more related products, see Axygen Product Selection Guide (CLS-A-PSG-001) at www.corning.com/lifesciences.



Australian Distributors:

FB AUSTRALIA

Free Call: 1800 066 077

Email: info@fisherbiotec.com

Web: fisherbiotec.com.au

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences
www.corning.com/lifesciences

NORTH AMERICA

t 800.492.1110
t 978.442.2200

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CSEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400