


Calibration kit, d18mm

DESCRIPTION

McFarland Turbidity Standards, Ø18mm

SPECIFICATIONS

Diameter	18 mm
----------	-------



CAT. NUMBER

70900

Ø18mm

ACCESSORIES



DEN-1

BS-050102-AAF

Densitometer (suspension turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ($0 - 180 \times 10^7$ cells/ml).

Densitometer provides the opportunity to measure solution ...

[read more](#)



DEN-1B

BS-050104-AAF

Densitometer (suspension turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ($0 - 180 \times 10^7$ cells/ml).

Densitometers provide the opportunity to measure solution turbidity in ...

[read more](#)