

Verification set for Abs

DESCRIPTION

Certified reference material, neutral density glass filter set of 4 Abs calibration points - 0.3532, 1.0512, 2.0425, 2.927 (the values may vary slightly from batch to batch)

CAT. NUMBER

BS-050109-AK

Neutral density glass filter set of 4

ACCESSORIES



DEN-600
BS-050109-AAA
Photometer

DEN-600 is a compact, portable, rechargeable battery powered photometer. DEN-600 comprises of 600 nm wavelength optical system, which enables to apply - 1) OD600 method that estimated total number of ...

[read more](#)