

# HTP-1, holder

## DESCRIPTION

HTP-1, holder for temperature probe



## CAT. NUMBER

BS-010309-FK

HTP-1, holder for temperature probe

## ACCESSORIES



**MMS-3000**  
BS-010305-AAF

High speed magnetic stirrer

**MMS-3000** are compact magnetic stirrers with stainless steel working surface. Units provide stirring of liquids with rotation speed of magnetic element up to 3000 rpm. Up to date it is ...

[read more](#)



**MSH-300**  
BS-010302-OAA  
Magnetic Stirrer with hot plate

**MSH-300** is magnetic stirrer of the new generation. Enclosures of stirrers are made of metal painted with powder enamel chemically resistant to acids and alkali. The stirrers are equipped with ...

[read more](#)



**Intelli-Stirrer MSH-300i**  
BS-010309-AAA  
Magnetic Stirrer with hot plate

**Intelli-Stirrer MSH-300i** is magnetic stirrer of the new generation. Enclosures of stirrers are made of metal painted with powder enamel chemically resistant to acids and alkali. The stirrers are equipped ...

[read more](#)