

MagSorb-16, Magnetic Rack for Manual Nucleic Acid Extraction

DESCRIPTION

Biosan presents the complete line of necessary instrumentation to utilize magnetic bead extraction kits and protocols from various manufacturers and meet the most demanding user requirements.

The one of the foundations of this line is **MagSorb-16** which is a magnetic rack that easily accommodates up to 16 single use tubes (1.5–2 ml).

Magnetic beads based NA extraction methods offer quick and efficient separation of genetic material from cellular leftovers. Wide variety of commercially available extraction kits manufactured by companies like ThermoFisher Scientific, Promega, Qiagen etc. are available in the market.

Our **MagSorb-16** magnetic rack can be used together with any manufacturers magnetic beads based NA extraction reagents, that allow working in 1.5- 2.0 ml tube format.



SPECIFICATIONS

Number of places in stand	16
Tube's volume	1.5 - 2 ml
Overall dimensions (W×D×H)	160x50x50 mm
Weight	0.24 kg

CAT. NUMBER

BS-010601

MagSorb-16, magnetic rack