

# Axygen® Axyspin Refrigerated Microcentrifuge



The Axygen Axyspin refrigerated microcentrifuge is a high performance refrigerated microcentrifuge designed to quickly process samples for applications such as nucleic acid and protein preparation, purifications, extractions, and pelleting.

- ▶ Powerful refrigeration system cools to 4°C in 8 minutes
- ▶ 24 x 1.5/2.0 mL aluminum rotor included
- ▶ Exceptionally quiet and compact

## Axygen Axyspin Refrigerated Microcentrifuge

The Axygen Axyspin refrigerated microcentrifuge combines a powerful yet silent, brushless drive with an efficient cooling system, all in a compact footprint. Brushless drive technology has been fine-tuned in the Axyspin refrigerated microcentrifuge to deliver fast acceleration and deceleration of the solid aluminum rotor. The uniquely designed 24-place rotor allows easy access to the tops of sample tubes. Starting from room temperature, the chamber of the Axyspin refrigerated microcentrifuge cools to 4°C in just 8 minutes. The control panel features a backlit LCD display, which allows users to easily view running parameters while the large control knob facilitates digital setting of rotor speed, run time, and temperature. These high performance features allow for faster sample processing and less down time.

The rotor of the Axyspin refrigerated microcentrifuge can accommodate 1.5 and 2.0 mL Axygen microcentrifuge tubes, including the most popular catalog numbers (MCT-150-C and MCT-200-C).

## Specifications

Speed Range	500 to 13,500 rpm
Maximum RCF	17,135 x g
Maximum Capacity	24 x 1.5/2.0 mL
Timer	0.5 to 99 minutes, or continuous 0.5 to 99 minutes, or continuous "Quick" button for momentary operation
Temperature Range	-10° to 40°C
Dimensions (W x D x H)	10.9 x 17.75 x 9.75 in. (27.7 x 45 x 24.75 cm)
Weight	44 lbs. (20 kg)
Electrical	100V-230V, 50/60 Hz

## Ordering Information

### Axygen Axyspin Refrigerated Microcentrifuge

Corning Cat. No.	Voltage	Plug Type	Qty/Cs
601-05-031	120V	US	1
601-05-021	230V	EU/UK	1
601-05-011	100V	US	1

### Axygen Axyspin MaxyClear SnapLock Microcentrifuge Tubes

Corning Cat. No.	Description	Qty/Pk	Qty/Cs
MCT-150	1.5 mL	500	5,000
MCT-200	2.0 mL	500	5,000



Australian distributors:  
 Fisher Biotec Australia  
 free call: 1800 066 077  
 email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
 web: [www.fisherbiotec.com](http://www.fisherbiotec.com)