





Reagent Reservoirs

Axygen® disposable reservoirs are designed to accommodate all multichannel pipettors. Four corners and a sharp V-bottom shape ensure maximum sample retrieval.

Disposable Reservoirs	Maximum Volume	Maximum Volume/Well	Color	Standard Cat. No.	Description
	25 mL	25 mL	White	RES-2CV-25	25 mL reservoir, divided
			White	RES-2CV-25-S*	25 mL reservoir, divided, sterile
			White	RES-2CV-25-SI	25 mL reservoir, divided, individually wrapped, sterile
Quantity: 100 reservoirs/case; * 5 reservoirs/pack, 40 packs/case					
	25 mL	25 mL	White	RES-V-25	25 mL reservoir
			White	RES-V-25-S*	25 mL reservoir, sterile
			White	RES-V-25-SI	25 mL reservoir, individually wrapped, sterile
Quantity: 100 reservoirs/case; * 5 reservoirs/pack, 40 packs/case					
	50 mL	50 mL	White	RES-V-50	50 mL reservoir
			White	RES-V-50-S*	50 mL reservoir, sterile
			White	RES-V-50-SI	50 mL reservoir, individually wrapped, sterile
Quantity: 100 reservoirs/case; * 5 reservoirs/pack, 40 packs/case					
	100 mL	100 mL	White	RES-V-100	100 mL reservoir
			White	RES-V-100-S*	100 mL reservoir, sterile
			White	RES-V-100-SI	100 mL reservoir, individually wrapped, sterile
Quantity: 100 reservoirs/case; * 5 reservoirs/pack, 40 packs/case					

