

Axygen® PlateMax® Heat Sealing Films

Axygen PlateMax microplate heat sealing films are designed for applications ranging from PCR and Real-time PCR to storage. They can be used in a wide range of temperatures making them suitable for nearly any application.

| Cat. No. | Material | Thickness (µm) | Dimensions (mm) | Min. Working Temp. | Max. Working Temp. | Application | Packaging |
|------------------------------|----------------------|----------------|-----------------|--------------------|--------------------|---|-----------------------|
| Paper-backed | | | | | | | |
| MF-300 | Litto/ethylene/vinyl | 85 | 127.2 x 79.5 | -80°C | 104°C | Storage | Bulk, 100/pk |
| UltraClear, Permanent | | | | | | | |
| HS-100-QPCR | Polyester | 52 | 125 x 78 | -80°C | 110°C | qPCR, fluorescence, colorimetric measurements | Bulk, 100/pk, 5 pk/cs |
| UltraClear, Peelable | | | | | | | |
| HS-150-QPCR | Polyester | 100.5 | 125 x 78 | -80°C | 80°C* | qPCR, fluorescence, colorimetric measurements | Bulk, 100/pk, 5 pk/cs |
| Peelable | | | | | | | |
| HS-200 | Polyester/aluminum | 70 | 125 x 78 | -80°C | 90°C* | Compound, storage, PCR | Bulk, 100/pk, 5 pk/cs |
| HS-400 | Polyester/aluminum | 70 | 125 x 78 | -80°C | 40°C | Compound storage | Bulk, 100/pk, 5 pk/cs |
| HS-300 | Aluminum | 20 | 125 x 78 | -80°C | 120°C | Storage | Bulk, 100/pk, 5 pk/cs |

*110°C with pressurized PCR-heated lids.




Axygen® PlateMax® Heat Sealing Films, Rolls

Axygen PlateMax microplate heat sealing films are designed for applications ranging from PCR and Real-time PCR to storage. They can be used in a wide range of temperatures making them suitable for nearly any application.



Pierceable

| Cat. No. | Material | Dimensions | Application | Compatibility | Packaging |
|--------------|----------|---------------|-------------|---------------|-----------------|
| HSF-AS-300-L | Aluminum | 115 mm x 500M | – | Agilent | Individual roll |
| HSF-AS-300-P | Aluminum | 78 mm x 610M | – | ALPS | Individual roll |



CyclerSeal, Peelable

| Cat. No. | Material | Dimensions | Application | Compatibility | Packaging |
|----------|------------------------|---------------|-------------|---------------|-----------------|
| HSF-TS-L | Polyester/ aluminum | 115 mm x 500M | PCR/Storage | Agilent | Individual roll |



Sealing Film Roller

Ideal for application of all Axygen PlateMax sealing films

| Cat. No. | Description | Packaging |
|---------------|---------------------|-----------|
| PCR-SP-ROLLER | Sealing film roller | 1/cs |