

Data Sheet

Biotin Elution Buffer BX (5 x Buffer BX)

Cat. No.: 2-1040-050

Lot No.: 1019-

Last date of revision

June 2014

Version 1040-1

Description	Biotin Elution Buffer BX comes with a high Biotin concentration for competitive displacement of purified recombinant <i>Strep</i> -tag [®] or Twin- <i>Strep</i> -tag [®] fusion proteins from immobilized Strep-Tactin XT . Buffer is recommended for efficient elution in combination with MagStrep "type3" XT beads . Please note that after elution with this buffer Strep-Tactin resin cannot be regenerated.
5 x concentrated solution	Buffer BX (Biotin elution buffer BX, 5x, 50 ml) 500 mM Tris-Cl pH 8.0 750 mM NaCl 5 mM EDTA 50 mM D-biotin
Working solution	Dilute 1 volume of concentrated solution (5x) with 4 volumes dest. water
Storage	4°C for concentrated and working solution
Stability	6 months after shipping
Shipment	RT

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.