



- Keep Containers & Glassware In Place
- Prevent Dangerous And Costly Spills
- Wipe Clean To Restore Adhesive
- Available In 2 Sizes
- Reusable Thousands of Times
- Ideal For Glassware Trolleys

### Keep Containers Where They Belong - On The Bench

Accidental spills in the lab can be dangerous and costly, but in many cases could be easily prevented with Ratek's unique washable adhesive work mats.

The work mats are adhesive on both sides and can be stuck indefinitely to any clean, dry and flat surface, yet are removable should they need to be relocated.



The unique high-tack adhesive has enough strength to prevent flasks and flat bottomed containers being accidentally tipped over, yet will still allow them to be removed when the time is right.

Once the mat becomes dirty, simply wipe clean with a damp cloth to restore the adhesive and allow it to dry, it's as simple as that.

Ideal for laboratory benches, glassware trolleys and sample preparation areas, the mats are available in 2 sizes to suit most applications. Mats can be easily trimmed to size where required.

### RELATED PRODUCTS

Product Code	Description
AT1	Adhesive mat tray to suit GM1 3D mixer
AT5	Adhesive mat tray to suit medium mixers
AT10	Adhesive mat tray to suit medium incubators

### ORDER CODES

Product Code	Description
AWM6X6	Washable adhesive work mat - 60 x 60mm
AWM33X16	Washable adhesive work mat - 332 x 166mm



AWM33X16



AWM6X6

