



SAFESEAL Microcentrifuge Tubes




Clearly the strongest tube, SafeSeal Microcentrifuge Tubes incorporate superior clarity with phenomenal strength.



- Superior clarity
- RCF rating of 24,000 x g
- Temperature rating from -80°C to +121°C
- Fit all standard rotors
- Autoclavable
- **WILL NOT POP OPEN DURING BOILING OR FREEZING**

Glass-like clarity allows easy viewing of sample and pellets. Ideal for high-speed centrifugation, proprietary homopolymer polypropylene tubes are rated up to 24,000 x g. Trust SafeSeal microcentrifuge tubes with your most valuable samples. A positive seal design allows repeated cap closures that will not pop during boiling. All tubes feature external graduation marks and frosted writing surfaces. Available in natural and five additional colors. Side gated tubes offer a clearer pellet view and increased strength over traditional style tubes.





SAFESEAL Microcentrifuge Tubes


0.65 ml SAFESEAL Microcentrifuge Tubes



Cat. No.	Color	Packaging	Quantity
11160	Natural	Non-Sterile, Bulk Pack (Box)	1,000 Tubes/Bag, 10 Packs/Case
11150	Natural	Non-Sterile, Bulk Bag	1,000 Tubes/Bag, 10 Bags/Case
11120	Red	Non-Sterile, Bulk Bag	1,000 Tubes/Bag, 10 Bags/Case
11060	Blue	Non-Sterile, Bulk Bag	1,000 Tubes/Bag, 10 Bags/Case
11080	Green	Non-Sterile, Bulk Bag	1,000 Tubes/Bag, 10 Bags/Case
11190	Yellow	Non-Sterile, Bulk Bag	1,000 Tubes/Bag, 10 Bags/Case
11170	Purple	Non-Sterile, Bulk Bag	1,000 Tubes/Bag, 10 Bags/Case
11210	Rainbow Pack*	Non-Sterile, Bulk Bag	1,000 Tubes/Bag, 10 Bags/Case

*Rainbow pack consists of 200 tubes each of red, blue, green, yellow and purple.


1.7 ml SAFESEAL Microcentrifuge Tubes



11500	Natural	Non-Sterile, Bulk Pack (Box)	500 Tubes/Bag, 10 Packs/Case
11510	Natural	Non-Sterile, Bulk Bag	500 Tubes/Bag, 10 Bags/Case
11490	Red	Non-Sterile, Bulk Bag	500 Tubes/Bag, 10 Bags/Case
11430	Blue	Non-Sterile, Bulk Bag	500 Tubes/Bag, 10 Bags/Case
11450	Green	Non-Sterile, Bulk Bag	500 Tubes/Bag, 10 Bags/Case
11560	Yellow	Non-Sterile, Bulk Bag	500 Tubes/Bag, 10 Bags/Case
11540	Purple	Non-Sterile, Bulk Bag	500 Tubes/Bag, 10 Bags/Case
11590	Rainbow Pack*	Non-Sterile, Bulk Bag	500 Tubes/Bag, 10 Bags/Case

*Rainbow pack consists of 100 tubes each of red, blue, green, yellow and purple.

2.0 ml SAFESEAL Microcentrifuge Tubes




12030	Natural	Non-Sterile, Bulk Pack (Box)	400 Tubes/Bag, 10 Packs/Case
12000	Natural	Non-Sterile, Bulk Bag	400 Tubes/Bag, 10 Bags/Case
11980	Red	Non-Sterile, Bulk Bag	400 Tubes/Bag, 10 Bags/Case
11920	Blue	Non-Sterile, Bulk Bag	400 Tubes/Bag, 10 Bags/Case
11940	Green	Non-Sterile, Bulk Bag	400 Tubes/Bag, 10 Bags/Case
12040	Yellow	Non-Sterile, Bulk Bag	400 Tubes/Bag, 10 Bags/Case
12020	Purple	Non-Sterile, Bulk Bag	400 Tubes/Bag, 10 Bags/Case
12060	Rainbow Pack*	Non-Sterile, Bulk Bag	400 Tubes/Bag, 10 Bags/Case

*Rainbow pack consists of 80 tubes each of red, blue, green, yellow and purple.



Low Binding Microcentrifuge Tubes





- RCF rating of 24,000 x g
- Temperature rating from -80°C to +121°C
- Low protein binding, ideal for Proteomics
- Non-reactive
- Will not pop during boiling or freezing
- Complete sample recovery
- Autoclavable

	Cat. No.	Color	Packaging	Quantity
 <p>0.65 ml Low Binding Microcentrifuge Tubes</p>	11300	 Natural	Non-Sterile, Pack Box	500 Tubes/Pack, 10 Packs/Case
 <p>1.7 ml Low Binding Microcentrifuge Tubes</p>	39640T	 Natural	Non-Sterile, Pack Box	250 Tubes/Pack, 10 Packs/Case
 <p>2.0 ml Low Binding Microcentrifuge Tubes</p>	12160	 Natural	Non-Sterile, Pack Box	200 Tubes/Pack, 10 Packs/Case

Dolphin Microcentrifuge Tubes




- Ideal design for pelletizing samples
- RCF rating of 15,000 x g
- Temperature rating from -80°C to +121°C
- Will not pop during boiling or freezing
- Fit all standard rotors
- Autoclavable

	Cat. No.	Color	Packaging	Quantity
 <p>2.0 ml Dolphin Microcentrifuge Tubes</p>	16140	 Natural	Non-Sterile, Bulk	1,000 Tubes/Bag, 5 Bags/Case












TabTop Microcentrifuge Tubes



Vertical tab eases handling of tubes in rotors, racks, heat blocks and baths. Hazardous aerosols are vented out the back of the tube and away from the researcher by pressing the levered TabTop. 0.5 ml tubes feature two sectors (15 μ l capacity) and 1.5 ml tubes feature four sectors (20 μ l capacity) for reactant loading prior to centrifugation.

- RCF rating of 22,000 x g
- Temperature rating from -80°C to +121°C
- Will not pop during boiling or freezing
- Autoclavable

	Cat. No.	Color	Packaging	Quantity
0.5 ml TabTop Microcentrifuge Tubes 	16890	 Natural	Non-Sterile, Bulk	1,000 Tubes/Bag, 10 Bags/Case
	16130	 Natural	Non-Sterile, Bulk	500 Tubes/Bag, 10 Bags/Case
1.5 ml TabTop Microcentrifuge Tubes 	16490	 Red	Non-Sterile, Bulk	500 Tubes/Bag, 10 Bags/Case
	16470	 Blue	Non-Sterile, Bulk	500 Tubes/Bag, 10 Bags/Case
	16480	 Green	Non-Sterile, Bulk	500 Tubes/Bag, 10 Bags/Case
	16510	 Yellow	Non-Sterile, Bulk	500 Tubes/Bag, 10 Bags/Case
	16500	 Purple	Non-Sterile, Bulk	500 Tubes/Bag, 10 Bags/Case

μ ltra Racks



- Supports up to 48 tubes in an alpha-numeric configuration
- Non-skid rubber feet hold the rack in place
- Manufactured from 100% recycled polypropylene

	Cat. No.	Fits Tube Size	Capacity	Quantity
μltra Rack Microcentrifuge Tube and Vial Racks	15170	0.5 ml, 0.65 ml	48 Tubes	6 Racks/Case
	15190	1.5 ml, 1.7 ml	48 Tubes	6 Racks/Case




TwistTop Vials



- Superior clarity.
- RCF rating of 17,000 x g
- Temperature rating from -80°C to + 121°C
- Two full turns of the cap create a reliable seal
- Autoclavable
- Fit all standard rotors
- Ethylene/propylene O-ring in cap for complete sealing
- Taller caps for easier handling



	Cat. No.	Cap Color	Packaging	Quantity
0.65 ml TwistTop Vial, Skirted 	15530	N/A (Caps Not Included)	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
	12330	 Natural	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
	12400	 Natural	Pre-Sterilized, Bulk	50 Vials/Bag, 10 Bags/Case
	12370	 Red	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
	12440	 Red	Pre-Sterilized, Bulk	50 Vials/Bag, 10 Bags/Case
	12350	 Blue	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
	12420	 Blue	Pre-Sterilized, Bulk	50 Vials/Bag, 10 Bags/Case
	12360	 Green	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
	12430	 Green	Pre-Sterilized, Bulk	50 Vials/Bag, 10 Bags/Case
	12380	 Yellow	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
12450	 Yellow	Pre-Sterilized, Bulk	50 Vials/Bag, 10 Bags/Case	
12340	 White	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case	
12410	 White	Pre-Sterilized, Bulk	50 Vials/Bag, 10 Bags/Case	
12390	 Rainbow Pack*	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case	

*200 caps of each color are bagged in each case separately.






TwistTop Vials



	Cat. No.	Cap Color	Packaging	Quantity	
1.7 ml TwistTop Vial, Skirted 	15570	N/A (Caps Not Included)	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case	
	12470 12540	Natural Natural	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
	12510 12580	Red Red	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
	12490 12560	Blue Blue	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
	12500 12570	Green Green	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
	12520 12590	Yellow Yellow	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
	12480 12550	White White	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
	12530	Rainbow Pack*	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case	
	1.7 ml TwistTop Vial, Conical 	15590	N/A (Caps Not Included)	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
		12610 12680	Natural Natural	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
		12650 12880	Red Red	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
		12630 12700	Blue Blue	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
12640 12710		Green Green	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
12660 12890		Yellow Yellow	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
12620 12690		White White	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case	
12670		Rainbow Pack*	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case	

*200 caps of each color are bagged in each case separately.






TwistTop Vials

	Cat. No.	Cap Color	Packaging	Quantity
2.0 ml TwistTop Vial, Skirted 	15610	N/A (Caps Not Included)	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
	12910 12980	Natural Natural	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
	12950 13020	Red Red	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
	12930 13000	Blue Blue	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
	12940 13010	Green Green	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
	12960 13030	Yellow Yellow	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
	12920 12990	White White	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
	12970	Rainbow Pack*	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case
	2.0 ml TwistTop Vial, Conical 	15630	N/A (Caps Not Included)	Non-Sterile, Bulk
13050 13120		Natural Natural	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
13090 13160		Red Red	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
13070 13140		Blue Blue	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
13080 13150		Green Green	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
13100 13170		Yellow Yellow	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
13060 13130		White White	Non-Sterile, Bulk Pre-Sterilized, Bulk	500 Vials/Bag, 2 Bags/Case 50 Vials/Bag, 10 Bags/Case
13110		Rainbow Pack*	Non-Sterile, Bulk	500 Vials/Bag, 2 Bags/Case

*200 caps of each color are bagged in each case separately.





TwistTop Vials

	Cat. No.	Cap Color	Packaging	Quantity
TwistTop Vial Caps 	15490	 Natural	Non-Sterile, Bulk	500 Caps/Bag, 2 Bags/Case
	15410	 Red	Non-Sterile, Bulk	500 Caps/Bag, 2 Bags/Case
	15470	 Blue	Non-Sterile, Bulk	500 Caps/Bag, 2 Bags/Case
	15460	 Green	Non-Sterile, Bulk	500 Caps/Bag, 2 Bags/Case
	15400	 Yellow	Non-Sterile, Bulk	500 Caps/Bag, 2 Bags/Case
	15480	 White	Non-Sterile, Bulk	500 Caps/Bag, 2 Bags/Case
	15510	 Rainbow Pack*	Non-Sterile, Bulk	200 Caps/Bag, 5 Bags/Case

*200 caps of each color are bagged in each case separately.



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com