

TrueLine



TrueLine

TrueLine offers a wide range of quality laboratory plastic products for use in many types of life science laboratories. Whether you need products for culturing cells, carrying out assays or doing PCR work, TrueLine offers the most popular sizes of the most used products in many labs.

TrueLine Cell Culture Dishes

- Clear polystyrene
- Surface for cell attachment
- Sterile
- Stackable
- Cover allows gas exchange



Cat. No.	Product Description	Qty/Pk	Qty/Pk
TR4000	10x35mm	20	500
TR4001	15x60mm	20	500
TR4002	20x100mm	20	500
TR4003	25x150mm	5	60



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com