

TrueLine



TrueLine

TrueLine offers a wide range of quality laboratory plastic products for use in many types of life science laboratories. Whether you need products for culturing cells, carrying out assays or doing PCR work, TrueLine offers the most popular sizes of the most used products in many labs.

TrueLine Centrifuge Tubes

- Polypropylene
- Sterile and non sterile options
- Bulk and styrofoam racked options
- Flat, bright green colored caps
- RNase-free and DNase-free
- Non pyrogenic



Cat. No.	Volume	Sterile	Max RCF	Packaging	Qty/Pk	Qty/Cs
TR2000	15ml	Yes	8,000	Rack	50	500
TR2001	15ml	Yes	8,000	Bulk	50	500
TR2002	15ml	No	8,000	Bulk	125	1250
TR2003	50ml	Yes	10,000	Rack	25	500
TR2004	50ml	Yes	10,000	Bulk	50	500
TR2005	50ml	No	10,000	Bulk	60	600



Australian distributors:
 Fisher Biotec Australia
 free call: 1800 066 077
 email: info@fisherbiotec.com
 web: www.fisherbiotec.com