

Compact Scales

HT Series

With 3 color bars



How would you use the three color bars?

At a school, you may want to use different colors for different classes or classrooms. In food service, you could mark scales that require special hygienic attention with a particular color. The HT Series are compact scales that challenge your creativity!

There is more...



A&D
A&D Australasia
"Excellence in Measurement"
www.andaustralasia.com.au



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com

Compact Scales **HT** Series

Large LCD Display

Clearly visible, 13.5 mm-high characters

Long Battery Life

Approx. 450 hours with alkaline batteries

Multiple Weighing Units

Selectable from g, oz, ozt, lb-oz*, tael, and tola

Auto Power Off Function

Automatically turns off the display after 5 minutes of inactivity

Tare Function

Offsets the weight of a container with just one press of a button

Security Slot

Slot in the rear for an anti-theft lock

*lb-oz not available for the HT-300



Storage Box

Penholders in Two Corners

Name Plate Holder



It is possible to perform weighing without taking out the scale.

The stackable case takes up less space for storage.



Specifications

Model	HT-300	HT-500	HT-3000	HT-5000
Weighing capacity	310 g	510 g	3100 g	5100 g
Minimum weighing value	0.1 g		1 g	
Power supply	4 AA batteries or an AC adapter (optional)			
Pan size	132 (W) × 130 (D) mm			
Dimensions	195 (W) × 136 (D) × 44 (H) mm			



Stainless Steel Pan (optional)



AC Adapter (optional)

Options Stainless steel pan (HT-10)
AC adapter

A&D
A&D Australasia

ACN 26 007 556 809

Head Office
32 Dew Street
THEBARTON South Australia 5031
Telephone (08) 8301 8100
Facsimile (08) 8352 7409

"Excellence in Measurement"

Victorian Office
39 Bakehouse Road
KENSINGTON Victoria 3031
Telephone (03) 9371 1555
Facsimile (03) 9372 1193

New South Wales Office
Unit 4, 14 Abbott Road
SEVEN HILLS New South Wales 2147
Telephone (02) 9674 5466
Facsimile (02) 9674 2544

sales@andweighing.com.au
www.andweighing.com.au

Your Dealer