

# Compact Scales

# HT Series

With 3 color bars



## How would you use the three color bars?

At a school, you may want to use different colors for different classes or classrooms. In food service, you could mark scales that require special hygienic attention with a particular color. The HT Series are compact scales that challenge your creativity!

*There is more...*



**A&D**  
A&D Australasia  
"Excellence in Measurement"  
[www.andaustralasia.com.au](http://www.andaustralasia.com.au)



Australian distributors:  
Fisher Biotec Australia  
free call: 1800 066 077  
email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
web: [www.fisherbiotec.com](http://www.fisherbiotec.com)



# Compact Scales **HT** Series

## Large LCD Display

Clearly visible, 13.5 mm-high characters

## Long Battery Life

Approx. 450 hours with alkaline batteries

## Multiple Weighing Units

Selectable from g, oz, ozt, lb-oz\*, tael, and tola

## Auto Power Off Function

Automatically turns off the display after 5 minutes of inactivity

## Tare Function

Offsets the weight of a container with just one press of a button

## Security Slot

Slot in the rear for an anti-theft lock

\*lb-oz not available for the HT-300



## Storage Box

Penholders in Two Corners

Name Plate Holder



It is possible to perform weighing without taking out the scale.

The stackable case takes up less space for storage.



## Specifications

Model	HT-300	HT-500	HT-3000	HT-5000
Weighing capacity	310 g	510 g	3100 g	5100 g
Minimum weighing value	0.1 g		1 g	
Power supply	4 AA batteries or an AC adapter (optional)			
Pan size	132 (W) × 130 (D) mm			
Dimensions	195 (W) × 136 (D) × 44 (H) mm			

**Options** Stainless steel pan (HT-10)  
AC adapter



Stainless Steel Pan (optional)



AC Adapter (optional)



Australian distributors:  
Fisher Biotec Australia  
free call: 1800 066 077  
email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
web: [www.fisherbiotec.com](http://www.fisherbiotec.com)