


Axygen®



A Corning Brand

Conical Centrifuge Tubes

Rated to 12,000 RCF



- Printed graduations and writing area for sample identification
- Screw Caps can be closed with one hand
- Sterile options are also available
- RNase-/DNase-free
- Nonpyrogenic



| Cat. No. | Volume (mL) | Description | Packaging |
|---------------|-------------|--------------------------------------|-------------------------------|
| SCT-15ML-500 | 15 | Clear | 50/pk, 10 pk/cs |
| SCT-15ML-R-S | 15 | Clear, sterile, SAL 10 ⁻⁵ | 50/rack, 10 racks/cs |
| SCT-15ML-25-S | 15 | Clear, sterile, SAL 10 ⁻⁵ | 25/pk, 20 pk/cs |
| SCT-50ML-500 | 50 | Clear | 50/pk, 10 pk/cs |
| SCT-50ML-R-S | 50 | Clear, sterile, SAL 10 ⁻⁵ | 25/rack, 10 racks/pk, 2 pk/cs |
| SCT-50ML-25-S | 50 | Clear, sterile, SAL 10 ⁻⁵ | 25/pk, 20 pk/cs |



www.fisherbiotec.com

For all your life science research requirements

Free Call: 1800 066 077 Tel: +61-8-9382 3511 Email: info@fisherbiotec.com

