

# Corning<sup>®</sup> IMDM (Iscove's Modification of DMEM), 1X

CORNING

## Product Facts

Complete Certificate of Analysis available for each production lot, along with full formulation information



Corning Media products are manufactured by Mediatech, Inc.

### Formulation (Choose your specification)

With Earle's Salts and L-glutamine, without Ribonucleosides, Deoxyribonucleosides, and Sodium Bicarbonate

[+]Earles Salts; [-]Ribonucleosides, Deoxyribonucleosides, L-glutamine

With Earles Salts; Ribonucleosides, Deoxyribonucleosides, and L-glutamine

### Size

50 mL, 500 mL, 1 L

### Shelf Life

12 months, 36 months

### Storage

2°C to 8°C

PLEASE **CONTACT US** FOR ORDERING INFORMATION  
OR RING FREECALL **1800 066 077**

**fb** Australian distributors:  
Fisher Biotec Australia  
free call: 1800 066 077  
email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
web: [www.fisherbiotec.com](http://www.fisherbiotec.com)