



What to look for in a BIOLOGICAL SAFETY CABINET

Key features for proven safety



Isocide™

Antimicrobial external surface paint impregnated with silver ions can inhibit microbial growth to reduce bioburden and improve sample safety.



DC ECM Blower

Direct Current Electronically Commutated Motor (DC ECM) produces less heat and less likely to be affected by power fluctuation. Users can have 70% more energy savings than with AC motor.



Dynamic Chamber™

Negative-pressure side walls, blower, and plenum help prevent contaminants from escaping.



Ultra Low Particulate Air (ULPA) Filter

Creates ISO 14644-1 Class 3 Work Zone with its high efficiency of >99.999% at 0.1 - 0.3 μm which is 10x more effective than HEPA Filter.

ESCO
WORLD CLASS. WORLDWIDE.



Airstream® Plus
Class II Biological Safety Cabinet



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com