



MORE THAN  
**20**  
YEARS  
IN BUSINESS  
Established 1997

# PROMOTIONAL OFFER



[www.fisherbiotec.com](http://www.fisherbiotec.com)

For all your life science requirements

free call **1800 066 077**

*Solutions through innovation*

## Purchase any Select Digital Dry Bath and get at least 15% discount off list pricing

Single, dual and four block capacity

Temperature range: Ambient +5° to 150°C

USB connectivity enables data traceability

Dual display



from **\$620** ex-GST



This offer valid until ~~28 December 2018~~ – now extended to 31 January 2019

The Select BioProducts Dry Baths provide a broad temperature range up to 150°C, which makes them useful for a variety of applications in molecular biology, histology, clinical, environmental and industrial laboratories. They have been designed to incorporate the best combination of features including block capacity, dual displays, and an external temperature probe that gives life science researchers the confidence and efficiency they need when working with precious samples.

### ORDERING INFORMATION

*Cat. No	Description	Normal Price	Special Price
WS-SBD-1301-230V	Select Heat Digital Dry Bath, single block capacity, temperature range ambient +5° to 100°C, includes Au power cords	\$731	<b>\$620</b>
WS-SBD-1302-230V	Select Heat Digital Dry Bath, dual block capacity, temperature range ambient +5° to 100°C, includes Au power cords	\$919	<b>\$780</b>
WS-SBD-1301-230V	Select Heat Digital Dry Bath, four capacity, temperature range ambient +5° to 100°C, includes Au power cords	\$1,326	<b>\$1,125</b>

\*Quote special Catalog No. when ordering

Above prices in AUD and exclude GST

**For more information on this offer – free call 1800 066 077**

Terms and Conditions

All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Fisher Biotec reserves the right to withdraw offers at any time.

# Blocks for Digital Block Heater

Purchase any Block for only **\$150 ea.**

This offer valid until 31 January 2019



## SPECIFICATIONS

**Construction** – anodised, high grade, non-porous aluminium

**Thermometer well** – available in all blocks except those used with plates and slides

**Compatibility** – will work in any manufacturer's dry bath with an opening of 3.04 x 3.60 inches.

Select  
BioProducts  
For Life Science Research

## ORDERING INFORMATION

*Cat. No	Description	Special Price
WS-D1101	Solid block for machining (no holes)	<b>\$150</b>
WS-D1102	Block, 48 x 0.2 ml PCR tubes or 6 x 0.2 ml strips	<b>\$150</b>
WS-D1102A	Block, 20 x 2.0 ml tubes	<b>\$150</b>
WS-D1105	Block, 24 x 0.5 mL tubes	<b>\$150</b>
WS-D1105A	Block, 24 x 1.5 ml tubes	<b>\$150</b>
WS-D1106	Block, 35 x 6 mm tubes	<b>\$150</b>
WS-D1110	Block, 20 x 10 mm tubes	<b>\$150</b>
WS-D1112	Block, 20 x 12 mm tubes	<b>\$150</b>
WS-D1113	Block, 20 x 13 mm tubes	<b>\$150</b>
WS-D1115-TALL	Block, 12 x 15 ml centrifuge tubes	<b>\$150</b>
WS-D1116	Block, 12 x 15 or 16 mm tubes	<b>\$150</b>
WS-D1117	Block, 12 x 17 mm tubes	<b>\$150</b>
WS-D1120	Block, 6 x 20 mm tubes	<b>\$150</b>
WS-D1125	Block, 6 x 25 mm tubes	<b>\$150</b>
WS-D1150-TALL	Block, 5 x 50 ml centrifuge tubes	<b>\$150</b>
WS-D1100-PS	PopStopper for microtube blocks, D1102, D1102A, D1105, D1105A	<b>\$150</b>
WS-D1105-RACK	Rack accessory for 20 x 1.5 ml tubes, compatible with D1105A, pack of 3	<b>\$150</b>

\*Quote special Catalog No. when ordering

Above prices in AUD and exclude GST

For more information on this offer – free call 1800 066 077

Terms and Conditions

All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Fisher Biotec reserves the right to withdraw offers at any time.