



Compound Library Storage – Mini Tube System (MTS)

Axygen® Mini Tube Systems are available in 0.65 mL single mini tube or 1.1 mL single mini tube. They are also available in 8-cap or 12-cap strip format. Mini Tube System racks are compatible with Beckman Coulter, Tecan, and other automation workstations.

- ▶ Rack features a one-way lid design and alpha-numeric grid
- ▶ Molded-in condensation rings on the lid minimize cross contamination



Cat. No.	Description	Tube Volume (mL)	Pack Type	Packaging
MTS-11-C	Mini tube, single	1.1	Bulk	960/bag, 5 bags/cs
MTS-11-C-R	Mini tube, single	1.1	Racked	96/rack, 10 racks/pk, 5 pk/cs
MTS-11-C-R-S	Mini tube, single, sterile	1.1	Racked	96/rack, 10 racks/pk, 5 pk/cs



MTS-11-8-C	Mini tube, 8-strip	1.1	Bulk	120/bag, 5 bags/cs
MTS-11-8-C-R	Mini tube, 8-strip	1.1	Racked	12/rack, 10 racks/pk, 5 pk/cs
MTS-11-8-C-R-S	Mini tube, 8-strip, sterile	1.1	Racked	12/rack, 10 racks/pk, 5 pk/cs



MTS-11-12-C	Mini tube, 12-strip	1.1	Bulk	80/bag, 5 bags/cs
MTS-11-12-C-R	Mini tube, 12-strip	1.1	Racked	8/rack, 10 racks/pk, 5 pk/cs
MTS-11-12-C-R-S	Mini tube, 12-strip, sterile	1.1	Racked	8/rack, 10 racks/pk, 5 pk/cs



MTS-06-C	Mini tube, single	0.65	Bulk	960/bag, 5 bags/cs
MTS-06-C-R	Mini tube, single	0.65	Racked	96/rack, 10 racks/pk, 5 pk/cs
MTS-06-C-R-S	Mini tube, single, sterile	0.65	Racked	96/rack, 10 racks/pk, 5 pk/cs



MTS-8CP-C	Mini tube, 8-cap strip	–	Bulk	120/pk, 5 pk/cs
MTS-8CP-C-S	Mini tube, 8-cap strip, sterile	–	Bulk	120/pk, 5 pk/cs



MTS-12CP-C	Mini tube, 12-cap trip	–	Bulk	80/pk, 5 pk/cs
MTS-12CP-C-S	Mini tube, 12-cap strip, sterile	–	Bulk	80/pk, 5 pk/cs

