



MAESTROGEN



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com



TK-02

iTrack cooling system

Specification

Dimensions(mm)	215 x 177 x 127 *
Weight(kg)	2.8
Power	DC 5.0V, 6.0A
Scan Light Source	IR Light-Emitting Diode
Scan Frequency(Hz)	200 - 500
Detection Time(sec)	0.5 - 3
Cooling capacity	< 4°C with 96-well metal rack adapter
Control System	Android tablet
Compatible Pipette	Single-, 8-, & 12-channel (fixed space)
Compatible Tip	0.1 - 1,000 µL
Compatible Microplate/Tube (PCR, ELISA, storage...)	96-Well Microplate (height < 23.5 mm) 384-Well Microplate (height < 11 mm) 8-tube / 12-tube strips

* the dimension exclude tablet