




## Gel Loading Tips



Sorenson is *THE* trusted manufacturer of gel loading tips. Sorenson's MultiFlex™, MultiFlex E and MiniFlex gel loading tips have become the industry standard and have been used globally in genetic sequencing studies.

### MultiFlex Tips (1-200 µl):


- Universal fit on research-grade pipettors
- 83 mm length for access to bottom of narrow vessels and gel wells
- 5 µl graduation for volume QuickCheck
- Available in bulk (round) and racked (round and flat)

### MiniFlex Tips (0.1-10 µl):

- Universal fit on ultra-micro pipettors
- 45 mm length for access to bottom of narrow vessels and gel wells.
- 2 µl graduation for volume QuickCheck
- Available racked (round and flat)

### MultiFlex E Tips (0.5-10 µl):

- Fit Eppendorf and ultra-micro pipettors
- 5 µl graduation for volume QuickCheck
- Microcapillary is ideal for viewing small sample volumes
- Available racked (round and flat)



White racks – Round (0.57 mm O.D.) tips  
 Yellow racks – Flat (0.37 mm O.D.) tips  
 Blue racks – Flat (0.17 mm O.D.) tips
















Australian distributors:  
 Fisher Biotec Australia  
 free call: 1800 066 077  
 email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
 web: [www.fisherbiotec.com](http://www.fisherbiotec.com)



# Gel Loading Tips



	Cat. No.	Vol. Range	Orifice Shape	Orifice O.D.*	Graduations	Packaging	Pack Type**	Quantity	
<b>10 µl MiniFlex</b> 	15100	0.1-10 µl	Round	0.57 mm	2 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case	
	15110	0.1-10 µl	Round	0.57 mm	2 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case	
	15120	0.1-10 µl	Flat	0.37 mm	2 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case	
	15130	0.1-10 µl	Flat	0.37 mm	2 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case	
	17350	0.1-10 µl	Flat	0.17 mm	2 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case	
	17360	0.1-10 µl	Flat	0.17 mm	2 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case	
	<b>10 µl MultiFlex E</b> 	13920	0.5-10 µl	Round	0.57 mm	5 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case
		13930	0.5-10 µl	Round	0.57 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
		13890	0.5-10 µl	Flat	0.37 mm	5 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case
13900		0.5-10 µl	Flat	0.37 mm	5 µl	Racked, Sterile	200 Tips/Pack, 4 Packs/Case		
17330		0.5-10 µl	Flat	0.17 mm	5 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case	
17340		0.5-10 µl	Flat	0.17 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case	
<b>200 µl MultiFlex</b> 		28480	1-200 µl	Round	0.57 mm	5 µl	Bulk, Non-Sterile		1000 Tips/Bag, 10 Bags/Case
		13790	1-200 µl	Round	0.57 mm	5 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case
		13810	1-200 µl	Round	0.57 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case
	13760	1-200 µl	Flat	0.37 mm	5 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case	
	13770	1-200 µl	Flat	0.37 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case	
	17310	1-200 µl	Flat	0.17 mm	5 µl	Racked, Non-Sterile		200 Tips/Pack, 4 Packs/Case	
	17320	1-200 µl	Flat	0.17 mm	5 µl	Racked, Sterile		200 Tips/Pack, 4 Packs/Case	

\*.37 mm O.D. tips are for use with 0.4 mm gels; .17 mm O.D. tips are for use with 0.2 mm gels

\*\*Gel Loading Tip Packaging: an individual pack is a tip rack of 200 tips



Australian distributors:  
 Fisher Biotec Australia  
 free call: 1800 066 077  
 email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
 web: [www.fisherbiotec.com](http://www.fisherbiotec.com)