




BandPeeper



MaestroGen offers BandPeeper real-time electrophoresis LED illuminator for DNA gel fluorescent signal observation in water solution.



- Time saving for new sample electrophoresis trial & error.
- Cost saving for repeated trial & error.
- Accurate electrophoresis duration control.
- Compatible with Mupid-size tank.
- Strong blue LED: used with common nucleic acid dyes.
- High sensitivity (~2 ng).
- Anti-fog design on amber filter.
- No decay for DNA gel (UV flashlamp decays DNA gel).
- Competitive price.



Order Information

Cat. No.	Window size	Wavelength	Sensitivity	Filter surface	Compatibility
RT-172	150 x 60 (mm)	LED 470 nm	up to 2 ng	Anti-Fog treatment	Mupid exU series



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com