


# UltraSlim<sup>®</sup> UV

## TRANSILLUMINATOR



UltraSlim UV transilluminator is a single-handedly take transilluminator with a high sensitivity performance to 1ng. A new-design UV blocking metal cover which is easily detachable for access to the work surface. UltraSlim UV transilluminator is an especially for electrophoresis mini-gel size and easy to observation or gel-cutting.

**ACCURATE**


**FLEXIBILITY**

**EASE**

**SENSITIVITY**

**SAFE**

## Specification



Cat No. UUV-01

Dimensions (mm) 247 (D) x 132 (W) x 48 (H)

View Size (mm) 150 (D) x 80 (W)

Wavelength (nm) 302 nm

Light Source UV lamp (T5-6W)

Uniformity +/- 5 % CV

Blocking Cover Stainless steel metal frame  
[Removable]

Filter UltraSafe UV Blocking  
(cover 99.99% UV transmittance)

Power DC12V / 1.5A

Weight (Kg) 1.4

