



# UV DARKROOMS



#### LARGE CAPACITY

Ideal for small or  
large sample.



#### UV ABSORBER SHIELD

Protect the user from  
the UV light.



#### COMBINATION OF UV SOURCES

UVA, UVB and UVC can be  
combined at your convenience.



#### EASY ACCESS TO DARKROOM

Ideal for mineralogy  
or TLC plate.

## About UV Darkrooms

---

Our UV darkrooms provide a large effective capacity and a high detection sensitivity. They could offer any combination of UV sources, simultaneously or not. The UV darkroom CN-6 is supplied without lamps and allows different lighting possibilities according to the user's choice.

## Darkrooms - 15W

Model	Wavelength (nm)	Power (W)	Type
CN-15	Multichannel darkroom without tubes, chose between the set below		
Set N°1 of 4 tubes	365 / 312	2 x 15 / 2 x 15	UV
Set N°2 of 4 tubes	365 / 254	2 x 15 / 2 x 15	UV
Set N°3 of 4 tubes	312 / 254	2 x 15 / 2 x 15	UV
Set N°4 of 4 tubes	365	4 x 15	UV
Set N°5 of 4 tubes	312	4 x 15	UV
Set N°6 of 4 tubes	254	4 x 15	UV
CN-6	Darkroom without lamp, the wavelengths can be chosen, between 254 nm or 365 nm - VL-6 Models		



White-Light Bulb For Normal  
Observation



Removable Lamps That Can Be  
Used For Hand-Held Applications