



UV LAMPS



UV MASTER TECHNOLOGY

Highly concentrated UV radiations.



EASY TO HANDLE

Lamp stand or holder to add versatility.



LONG LIFETIME

Anodized aluminum housing.



ONDULEX REFLECTOR

Reflect the maximum UV output.

About UV Lamps

Our large choice of filtered and unfiltered lamps range from 6 to 3x40 w-watt.

Our filtered lamps are ideal for the fluorescence detection. The BL-Black-Light lamps provide a large light spectrum with 365nm UV and visible light. The UV germicidal 254nm lamps help to prevent contamination in research laboratories.

Lamps 6W - Filtered

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-6.LC	365 / 254	145 x 48	2 x 6	UV Fluorescence
VL-6.L	365	145 x 48	1 x 6	UV Fluorescence
VL-6.C	254	145 x 48	1 x 6	UV Fluorescence

Lamps 6W - Unfiltered

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-206.BL	365	No Filter	2 x 6	Black Light
VL-206.G	254	No Filter	2 x 6	Germicidal

Lamps 15W - Filtered

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-215.L	365	295 x 76	2 x 15	UV Fluorescence
VL-215-M	312	295 x 76	2 x 15	UV Fluorescence
VL-215-C	254	295 x 76	2 x 15	UV Fluorescence
VL-115.L	365	295 x 66	1 x 15	UV Fluorescence
VL-115-M	312	295 x 66	1 x 15	UV Fluorescence
VL-115-C	254	295 x 66	1 x 15	UV Fluorescence

Lamps 15W - Unfiltered

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-315-BL	365	No Filter	3 x 15	Black Light
VL-315-G	254	No Filter	3 x 15	Germicidal
VL-215-BL	365	No Filter	2 x 15	Black Light
VL-215-G	254	No Filter	2 x 15	Germicidal
VL-115-BL	365	No Filter	1 x 15	Black Light
VL-115-G	254	No Filter	1 x 15	Germicidal

Lamps 30W -Unfiltered

Model	Wavelength (nm)	Filter Size (mm)	Power (W)	Type
VL-330.BL	365	No Filter	3 x 30	Black Light
VL-330.G	254	No Filter	3 x 30	Germicidal
VL-230.BL	365	No Filter	2 x 30	Black Light
VL-230.G	254	No Filter	2 x 30	Germicidal
VL-130.BL	365	No Filter	1 x 30	Black Light
VL-130.G	254	No Filter	1 x 30	Germicidal

Irradiation Surface For A Height Of About 2m

Power (W)	1 Tube	2 Tubes	3 Tubes
15	1.5 x 2m	1.5 x 3m	1.5 x 4m
20	2m x 2m	2m x 3m	2m x 4m
30	2.5m x 2.5m	2.5m x 3.5m	2.5m x 4.5m
40	3m x 3m	3m x 4m	3m x 5m

Lamps Stands & Handlers

Models	Type	Compatibility
SVL-30	Stand	For VL-115 and VL-215
SVL-6	Stand	For VL-6 and VL-8
SMA	Handle	For VL-6 and VL-8
SMU	Bracket	For VL-6 and VL-8



Lamp Stand Or Holder To Add
Versatility



Highly Concentrated UV
Radiations



Reflect The Maximum
UV Output