



# UV TABLES



## UV MASTER TECHNOLOGY

Highly concentrated UV  
radiations.



## LONG LIFETIME FILTER

Unlimited life expectancy for  
312 & 365nm.



## STAINLESS STEEL FRAME

Long lifetime and  
no rust.



## ADJUSTABLE DUAL INTENSITY

Ideal to switch from short gel  
visualization to longer preparative work.

## About UV Tables

---

We provide first-class UV transilluminators for documentation and preparation in different dimensions and configurations. Our innovative Super-Bright filter stops all the visible light emitted by the tubes. As the UV tubes are no more visible, the Super-Bright improves dramatically the quality of gel visualization and documentation.

## Transilluminators 8W

| Model           | Wavelength (nm) | Filter Size (mm) | Power (W) | Type           | Intensity (uW/cm <sup>2</sup> ) |
|-----------------|-----------------|------------------|-----------|----------------|---------------------------------|
| ECX-F26.MX V1   | 312             | 210 x 260        | 5 x 8     | UV SuperBright | UV SuperBright                  |
| ECX-F26.M V1    | 312             | 210 x 260        | 5 x 8     | UV             | UV                              |
| ECX-F26.C V1    | 254             | 210 x 260        | 5 x 8     | UV             | UV                              |
| ECX-F20.L V1    | 365             | 200 x 200        | 5 x 8     | UV             | UV                              |
| ECX-F20.M V1    | 312             | 200 x 200        | 5 x 8     | UV             | UV                              |
| ECX-F20.C V1    | 254             | 200 x 200        | 5 x 8     | UV             | UV                              |
| ECX-F20.BLUE V1 | -               | 200 x 200        | -         | Blue Light     | Blue Light                      |
| ECX-F20.W V1    | -               | No Filter        | 5 x 8     | White Light    | White Light                     |
| ECX-F15.M       | 312             | 150 x 150        | 4 x 8     | UV             | UV                              |

## Transilluminators 8W - Multiband

| Model         | Wavelength (nm) | Filter Size (mm) | Power (W)     | Type           |
|---------------|-----------------|------------------|---------------|----------------|
| TCP-26.LMX V1 | 365 / 312       | 210 x 260        | 4 x 8 / 5 x 8 | UV SuperBright |
| TCP-26.LC V1  | 365 / 254       | 210 x 260        | 4 x 8 / 5 x 8 | UV             |
| TCP-26.MC V1  | 312 / 254       | 210 x 260        | 4 x 8 / 5 x 8 | UV             |

## Conversion Screen

| Model   | From This Wavelength (nm) | To This Wavelength (nm) | Type        |
|---------|---------------------------|-------------------------|-------------|
| FC-26WL | 312                       | -                       | White Light |