

## Rocker 35 EZ

The Rocker 35 EZ is the next generation large capacity rocker from Labnet. This upgraded Rocker 35 is designed for durability in laboratories requiring a large platform area and precise control over rocking parameters. It is backed by a 5 year warranty.

Both speed and time are easily selectable and the analog control has been upgraded to display discreet values. The tilt is now easily adjustable with knob on side. This system can handle up to 11 pounds or 5 kg of load so a double platform can be fully loaded without any variation to speed.

The unit comes in both single and dual platform models and includes a ribbed no-slip rubber mat to keep containers in place. This unit is also compatible with our sticky mat (I-5240) for applications requiring greater sample adhesion to the platform.



### SPECIFICATIONS

Speed range	10 - 120 rpm
Timer	10 - 120 min or continuous
Tilt range	Adjustable, $\pm 20^\circ$
Maximum Load	11 lb/5 kg
Ambient operating range	-10° to 70°C
Dimensions (W x D x H)	
Single platform	14.2 x 14.2 x 6.7 in/36 x 36 x 17 cm
Double platform	14.2 x 14.2 x 10.2 in/36 x 36 x 26 cm
Weight	
Single platform	14.3 lb/6.5 kg
Double platform	17.2 lb/7.8 kg
Electrical	120V~, 50/60 Hz or 230~, 50/60 Hz

CAT NO.	DESCRIPTION
S2035-E*	Rocker 35 EZ, single platform
S2035-DE*	Rocker 35 EZ, double platform

\*To order 230V units add -230V to the end of the catalog number. Includes EU and UK cords