

Labpette FX™ Fixed Volume Pipettes

Labpette FX pipettes are ideal for applications where specific volumes are used repeatedly and sterilization by autoclaving is required. Thirty-two pipettes cover the range from 2 µl to 1,000 µl with exceptional accuracy and reproducibility.

An ergonomically designed body and finger rest make these pipettes comfortable to hold, even during extended pipetting sessions. The smooth plunger action and unique "twist-eject" tip removal system require minimal force and reduce the risk of repetitive strain injuries.

Manufactured from lightweight, high impact plastic, the Labpette FX pipettes are extremely durable, offer excellent chemical resistance and retain their high level of accuracy and precision during repeated use. The Labpette FX can be autoclaved at 121°C.

Labpette FX pipettes are color coded for easy identification and proper tip selection. Each pipette is individually tested, supplied with a factory calibration certificate and covered by a limited lifetime warranty.



CAT NO.	VOLUME	ACCURACY	PRECISION	NON-FILTERED TIP
P3900-2	2 µl	±3.0	±1.5	200 µl
P3900-5	5 µl	±2.0	±1.0	200 µl
P3900-7	7 µl	±1.5	±0.7	200 µl
P3900-10	10 µl	±1.2	±0.5	200 µl
P3900-15	15 µl	±1.0	±0.5	200 µl
P3900-20	20 µl	±0.9	±0.4	200 µl
P3900-25	25 µl	±0.9	±0.4	200 µl
P3900-30	30 µl	±0.9	±0.3	200 µl
P3900-40	40 µl	±0.9	±0.3	200 µl
P3900-50	50 µl	±0.9	±0.3	200 µl
P3900-60	60 µl	±0.9	±0.3	200 µl
P3900-70	70 µl	±0.9	±0.3	200 µl
P3900-75	75 µl	±0.9	±0.3	200 µl
P3900-80	80 µl	±0.9	±0.3	200 µl
P3900-90	90 µl	±0.9	±0.3	200 µl
P3900-100	100 µl	±0.9	±0.3	200 µl
P3900-120	120 µl	±0.8	±0.3	200 µl
P3900-150	150 µl	±0.7	±0.3	200 µl
P3900-200	200 µl	±0.6	±0.2	200 µl
P3900-200A	200 µl	±0.6	±0.2	1,000 µl
P3900-220	220 µl	±0.6	±0.2	1,000 µl
P3900-250	250 µl	±0.6	±0.2	1,000 µl
P3900-300	300 µl	±0.6	±0.2	1,000 µl
P3900-400	400 µl	±0.6	±0.2	1,000 µl
P3900-450	450 µl	±0.6	±0.2	1,000 µl
P3900-500	500 µl	±0.6	±0.2	1,000 µl
P3900-600	600 µl	±0.6	±0.2	1,000 µl
P3900-700	700 µl	±0.6	±0.2	1,000 µl
P3900-750	750 µl	±0.6	±0.2	1,000 µl
P3900-800	800 µl	±0.6	±0.2	1,000 µl
P3900-900	900 µl	±0.6	±0.2	1,000 µl
P3900-1000	1,000 µl	±0.6	±0.2	1,000 µl
P3985-P	Universal carousel stand for 6 pipettes			



See page 34 for recommended tips

Other volumes available, see your Labnet Representative for details.