

Mini Water Bath

The Mini Water Bath has a 6 liter capacity to accommodate test tube racks or up to three 250 mL flasks. A powerful 550W heater allows for quick temperature adjustment. A hydraulic thermostat ensures accurate temperature control. The seamless, stainless steel chamber is corrosion resistant and easy to clean. A stainless steel gable cover helps to retain heat in the bath and control evaporation. A clip is included for suspending a thermometer in the bath.



SPECIFICATIONS	Mini Water Bath
Temperature range	Ambient +5° to 100°C
Capacity	6 liters
Controls	Analog
Exterior dimensions (W x D x H)	14 x 11 x 13 in 36 x 28 x 32 cm
Interior dimensions	32 x 17 x 18 cm 12.6 x 6.7 x 6.9 in
Weight	7.7 kg/17 lb
Electrical	120V~or 230~, 50/60 Hz

CAT NO.	DESCRIPTION
W1106A*	Mini Water Bath with stainless steel gable cover and thermometer clip 120V
W0020-100C	Non-toxic, PTFE coated thermometer fro -20° to 110°C, with 1° resolution, 305mm long, 76mm immersion

*To order 230V units add -230V to the end of the catalog number. 230V units include EU and UK cords