

HighPrep™ Viral DNA/RNA Kit

Rapid and reliable isolation of total viral nucleic acids from whole blood, serum, plasma, saliva and other body fluids



Description

MagBio's HighPrep™ Viral DNA/RNA Kit is designed for rapid and reliable isolation of total nucleic acid from whole blood, serum, plasma, saliva and other bodily fluids. High quality DNA is suitable for direct use in most downstream applications such as amplifications and enzymatic reactions. The kit can be adapted to most major liquid handling workstations in the market.



HighPrep™ Viral DNA/RNA Kit

Applications

Viral DNA and RNA isolation for:

- RT-qPCR; RT-PCR
- Virus detection, genotyping
- Viral load monitoring

[Click to view Protocol >>](#)

Benefits

- Rapid and reliable purification of viral nucleic acids
- Adaptable to various automated liquid handling workstations
- No toxic organic solvents

Ordering Information

Cat No.	Product	Preps
HPV-DR96	HighPrep™ Viral DNA/RNA Kit (96 preps)	96
HPV-DR96x4	HighPrep™ Viral DNA/RNA Kit (384 preps)	384



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com