



High Capacity and Extremely Versatile

The RM5 is a versatile heavy duty roller mixer capable of turning large drums, barrels or a carrier carousel loaded with large quantities of smaller containers.

Up to 20 kg balanced load, 4-30 RPM, digital control

By altering the fixing points in the carousel, different mixing actions can be obtained, including end-over-end, rolling as well as tilt and roll. Carousels to secure virtually any tube, flask, bag or tube rack can be supplied. The long-life polyurethane rollers are mounted on heavy duty shafts with precision ball bearings and are driven by a heavy duty geared brushless induction motor controlled by Ratek's precision drive system. The inbuilt electronic timer can be set for 0-99 minutes or 0-99 hours, memory buttons allow for up to 3 speed / time patterns to be stored via the large user-friendly control panel.

Ratek can design and manufacture a carousel wheel to suit your application and extend the RM5 to virtually any length as required. Ideal for mixing large quantities of laboratory flasks, paint cans, bottles, soil samples or test tubes, the RM5 is designed to work day-in day-out under heavy loads. The RM5 comes with built-in Ratek reliability you can expect from Australia's leading laboratory mixer brand.