

# Axygen® MultiRack Tip System

**AXYGEN**  
A Corning Brand

The Axxygen MultiRack tip system has been engineered to fit most of the pipettors in the market. It is an interchangeable, complete system offered in bulk, racked, filter, and reload packaging styles to meet the variety of needs in the lab setting. The MultiRack tip system offers standard and extended-length tips allowing volumes from 0.1  $\mu\text{L}$  to 1250  $\mu\text{L}$ . For racked and reload packaging styles, the tips are pre-loaded into easy-load, color-coded deck inserts. With a small benchtop footprint and reusable rack, this system minimizes environmental impact and lab space. Axxygen MultiRack tip system products are made in the USA and meet ISO standards. All tips are free of detectable RNase, DNase, DNA, and pyrogens.

## Features

- ▶ Reusable rack with reload systems
- ▶ Small footprint and stackable design for maximum bench space usage
- ▶ Extended tips offered in 10, 200, 300, and 1000  $\mu\text{L}$  sizes
- ▶ Color-coded racks for easy volume identification







**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

**CORNING**

**Corning Incorporated**  
*Life Sciences*