

Aerosol Barrier Tips

- ▶ Technologically advanced filter material allows free passage of air while blocking aerosol contamination.
- ▶ Filter membrane contains no leachable additives.
- ▶ Axxygen® Maxymum Recovery® surface improves accuracy and ensures maximum sample recovery.
- ▶ RNase-/DNase-free
- ▶ Nonpyrogenic
- ▶ Human gDNA- and PCR-inhibition-free

10 µL, Microvolume

Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Tip Rack	Pack Type	Packaging
TF-300	–	0.5 - 10	A	Bulk	1,000/pk, 10 pk/cs
TF-300-R-S	TF-300-L-R-S	0.5 - 10	A	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

10 µL, Microvolume, Extended

TXLF-10	–	0.5 - 10	B	Bulk	1,000/pk, 10 pk/cs
TXLF-10-R-S	TXLF-10-L-R-S	0.5 - 10	B	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

10 µL, Ultra Micro

TF-400	–	0.5 - 20	C	Bulk	1,000/pk, 10 pk/cs
TF-400-R-S	TF-400-L-R-S	0.5 - 20	C	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

20 µL, Ultra Micro

TF-420-R-S	TF-420-L-R-S	0.5 - 20	D	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs
------------	--------------	----------	---	-----------------	-------------------------------

20 µL, Universal

TF-20	–	0.5 - 20	D	Bulk	1,000/pk, 10 pk/cs
TF-20-R-S	TF-20-L-R-S	0.5 - 20	D	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

50 µL, Universal

TF-50	–	0.5 - 50	D	Bulk	1,000/pk, 10 pk/cs
TF-50-R-S	TF-50-L-R-S	0.5 - 50	D	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

100 µL, Universal

TF-100	–	1 - 100	D	Bulk	1,000/pk, 10 pk/cs
TF-100-R-S	TF-100-L-R-S	1 - 100	D	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

150 µL, Universal

TF-150	–	1 - 150	D	Bulk	1,000/pk, 10 pk/cs
TF-150-R-S	–	1 - 150	D	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

200 µL, Universal


TF-200	–	1 - 200	D	Bulk	1,000/pk, 10 pk/cs
TF-200-R-S	TF-200-L-R-S	1 - 200	D	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

200 µL, Wide Bore*

TF-205-WB-R-S	TF-205-WB-L-R-S	1 - 200	D	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs
---------------	-----------------	---------	---	-----------------	-------------------------------

*Wide bore tips are ideal for pipetting viscous materials, genomic DNA, or delicate cell cultures.

Tip Racks



Aerosol Barrier Tips (Continued)

300 µL, universal



Standard Surface Cat. No.	Maxyimum Recovery® Surface Cat. No.	Volume (µL)	Tip Rack	Rack Type	Packaging
TF-350	–	1 - 300	A	Bulk	1,000/pk, 10 pk/cs
TF-350-R-S	TF-350-L-R-S	1 - 300	A	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

1,000 µL, Universal



TF-1000	–	100 - 1,000	B	Bulk	1,000/pk, 5 pk/cs
TF-1000-R-S	TF-1000-L-R-S	100 - 1,000	B	Racked, sterile	100/rack, 10 racks/pk, 5 pk/cs

1,000 µL, Wide Bore



TF-1005-WB-R-S	TF-1005-WB-L-R-S	100 - 1,000	B	Racked, sterile	100/rack, 10 racks/pk, 5 pk/cs
----------------	------------------	-------------	---	-----------------	--------------------------------

1,100 µL, Beveled



TF-1200-C-R-S	TF-1200-C-L-R-S	100 - 1,100	C	Racked, sterile	96/rack, 12 racks/pk, 3 pk/cs
---------------	-----------------	-------------	---	-----------------	-------------------------------

1,250 µL, Matrix-style



MTXF-1250-C-R-S	–	15 - 1,250	D	Racked, sterile	96/rack, 12 racks/pk, 3 pk/cs
-----------------	---	------------	---	-----------------	-------------------------------

Tip Racks

