

Microcentrifuge Tubes and Screw Caps

Axygen® microcentrifuge tubes feature a screw cap that provides a tight secure seal. Screw cap tubes are available in a variety of colors with a conical or self-standing bottom. Sterile options are also available.

- ▶ Maximum G-force is 20,000 RCF.
- ▶ RNase-/DNase-free
- ▶ Nonpyrogenic

Tubes, with Cap, Sterile

All sterile screw cap microcentrifuge tubes are pre-assembled with non-looped caps with O-rings.



Cat. No.	Volume (mL)	Style	Color	Packaging
SCT-45-C-S	0.5	Conical, elongated, sterile	Clear	100/bag, 5 bags/pk, 8 pk/cs
SCT-050-C-S	0.5	Conical, sterile	Clear	100/bag, 5 bags/pk, 8 pk/cs
SCT-050-SS-C-S	0.5	Self-standing, sterile	Clear	100/bag, 5 bags/pk, 8 pk/cs
SCT-150-C-S	1.5	Conical, sterile	Clear	100/bag, 5 bags/pk, 8 pk/cs
SCT-150-SS-C-S	1.5	Self-standing, sterile	Clear	100/bag, 5 bags/pk, 8 pk/cs
SCT-200-C-S	2.0	Conical, sterile	Clear	100/bag, 5 bags/pk, 8 pk/cs
SCT-200-SS-C-S	2.0	Self-standing, sterile	Clear	100/bag, 5 bags/pk, 8 pk/cs

Colored-cap versions are also available in sterile format.

Tubes Only



Cat. No.	Volume (mL)	Style	Color	Packaging
ST-45	0.5	Conical	Clear	500/pk, 8 pk/cs
ST-45-SS	0.5	Self-standing	Clear	500/pk, 8 pk/cs
ST-050	0.5	Conical	Clear	500/pk, 8 pk/cs
ST-050-SS	0.5	Self-standing	Clear	500/pk, 8 pk/cs
ST-150	1.5	Conical	Clear	500/pk, 8 pk/cs
ST-150-SS	1.5	Self-standing	Clear	500/pk, 8 pk/cs
ST-200	2.0	Conical	Clear	500/pk, 8 pk/cs
ST-200-SS	2.0	Self-standing	Clear	500/pk, 8 pk/cs
ST-050-X	0.5	Conical	Amber	500/pk, 8 pk/cs
ST-050-SS-X	0.5	Self-standing	Amber	500/pk, 8 pk/cs
ST-060-SS	0.5	Self-standing	Clear	500/pk, 8 pk/cs
ST-150-X	1.5	Conical	Amber	500/pk, 8 pk/cs
ST-150-SS-X	1.5	Self-standing	Amber	500/pk, 8 pk/cs
ST-200-X	2.0	Conical	Amber	500/pk, 8 pk/cs
ST-200-SS-X	2.0	Self-standing	Amber	500/pk, 8 pk/cs

Microcentrifuge Tubes and Screw Caps (Continued)

Screw Caps, Non-looped, with O-ring



Cat. No.	Cap Color	Packaging
SCO-C	Clear	500/pk, 8 pk/cs
SCO-R	Red	500/pk, 8 pk/cs
SCO-O	Orange	500/pk, 8 pk/cs
SCO-Y	Yellow	500/pk, 8 pk/cs
SCO-G	Green	500/pk, 8 pk/cs
SCO-B	Blue	500/pk, 8 pk/cs
SCO-BR	Brown	500/pk, 8 pk/cs

Screw Caps, Looped, with O-ring



Cat. No.	Cap Color	Packaging
SCO-LP-C	Clear	500/pk, 8 pk/cs
SCO-LP-R	Red	500/pk, 8 pk/cs
SCO-LP-O	Orange	500/pk, 8 pk/cs
SCO-LP-Y	Yellow	500/pk, 8 pk/cs
SCO-LP-G	Green	500/pk, 8 pk/cs
SCO-LP-B	Blue	500/pk, 8 pk/cs
SCO-LP-BR	Brown	500/pk, 8 pk/cs