



Microvolume Tips

Microvolume tips are available with Axygen® Maxymum Recovery® surface technology. Maxymum Recovery improves recovery of costly samples and offers greater precision and accuracy.


► See tip compatibility tables, pages 83-88.

10 µL, Microvolume

Standard Surface Cat. No.	Maxymum Recovery Surface Cat. No.	Volume (µL)	Tip Rack	Pack Type	Packaging
T-300	T-300-L	0.5 - 10	A	Bulk	1,000/pk, 20 pk/cs
T-300-R	T-300-L-R	0.5 - 10	A	Racked	96/rack, 10 racks/pk, 5 pk/cs
T-300-R-S	T-300-L-R-S	0.5 - 10	A	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs
T-300-STK	T-300-L-STK	0.5 - 10	A	Stack rack	96/rack, 5 racks/stack, 10 stacks/cs

10 µL, Microvolume, Extended

TXL-10	TXL-10-L	0.5 - 20	B	Bulk	1,000/pk, 20 pk/cs
TXL-10-R	TXL-10-L-R	0.5 - 20	B	Racked	96/rack, 10 racks/pk, 5 pk/cs
TXL-10-R-S	TXL-10-L-R-S	0.5 - 20	B	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs



Close up of tip end



20 µL, Ultra Micro

T-400	T-400-L	0.5 - 20	C	Bulk	1,000/pk, 20 pk/cs
T-400-R	T-400-L-R	0.5 - 20	C	Racked	96/rack, 10 racks/pk, 5 pk/cs
T-400-R-S	T-400-L-R-S	0.5 - 20	C	Racked, sterile	96/rack, 10 racks/pk, 5 pk/cs

Tip Racks

