

# COVID-19 VIRUS DETECTION



## BENEFITS

- Safer handling and fewer bio-hazard precautions
- Viral RNA stable at room temperature for up to 21 days
- Any swab can be used
- Validated for use with most extraction kits

[LEARN MORE >>](#)

## iSWAB™ Microbiome Collection Kits



VALIDATED BY THE CDC FOR COLLECTION OF CORONAVIRUS SAMPLES, DUE TO ITS PROPERTY OF INACTIVATING THE VIRUS AND NOT REQUIRING COLD CHAIN SHIPMENT

Catalog No	Product Name	Volume
ISWAB-MB-1200	iSWAB-Microbiome Collection Kit, 1.0mL	1.0mL
ISWAB-MBF-1200	iSWAB-Microbiome Collection Kit with FecesCatcher, 1.0mL	1.0mL
ISM-T-1200-R	iSWAB-Microbiome Collection Tube rack, 1.0mL	1.0mL



## HighPrep™ Viral DNA/RNA Kits

RAPID AND RELIABLE MAGNETIC  
BEAD BASED ISOLATION OF VIRAL  
NUCLEIC ACIDS

## APPLICATIONS

Viral DNA and RNA isolation for:

- RT-QPCR, RT-PCR
- Virus detection, genotyping
- Viral load monitoring

[LEARN MORE >>](#)

## BENEFITS

- Rapid and reliable purification of viral nucleic acids
- Adaptable to various automated liquid handlers
- No toxic organic solvents



HighPrep™ Viral DNA/RNA Kit (HPV-DR-96)

Catalog No	Product Name
HPV-DR96	HighPrep™ Viral DNA/RNA Kit (96 Preps)
HPV-DR96x4	HighPrep™ Viral DNA/RNA Kit (384 Preps)

CONTACT US FOR PRICING - 1800 066 077

# COVID-19 VIRUS DETECTION

**AXYGEN.**

A Corning Brand

Your source of high quality plates and sealing films for RT-PCR machines

Brand	Model	Compatible Plates	Product Description	Case Qty
Applied Biosystems	QuantStudio™ qPCR Systems ViiA 7™ Real-Time PCR System	PCR-96-AB-C	200 µl, 96-well, elevated skirt	100/pk, 5pk/cs
		PCR-96M2-HS-C	200 µl, 96-well, half skirt, single notch corner	100/pk, 5pk/cs
		PCR-96-LP-AB-C (fast reaction)	100 µl, 96-well PCR microplate, low profile	5 sleeves/pk, 5pk/cs
		PCR-384M2-C	55 µl, 384-well, skirted	10/pk, 5pk/cs
Roche	LightCycler® 480 Real-Time PCR System	PCR-96-LC480-W-NF	100 µl, 96-well PCR microplate, white, no sealing film	10/pk, 5pk/cs
		PCR-384-LC480-W-NF	100 µl, 96-well PCR microplate, white, no sealing film	10/pk, 5pk/cs
All brands	All models	UC-500	70 µm Ultra Clear Pressure Sensitive Sealing Film for Real Time PCR, Nonsterile	100/pk, 5pk/cs


## Axygen® PCR Microplates

[LEARN MORE >>](#)



## Axygen® Sealing Films

[LEARN MORE >>](#)



## High quality Liquid Handling Automation Tips

Axygen® produces a range of high quality tips designed to function at similar or higher levels than the original manufacturer's tips. Tips are available for all major automation platforms including, Agilent, Beckman, Hamilton, Perkin Elmer and Tecan.

[LEARN MORE >>](#)



## Deep Well Plates and Reservoirs for Automation Solutions

Axygen® assay and deep well plates are used for sample collection, *in vitro* growth chamber, and long-term storage. Made from polypropylene, they offer excellent chemical resistance and will withstand temperatures down to -80°C.

[LEARN MORE >>](#)

