



MORE THAN  
**20**  
YEARS  
IN BUSINESS  
Established 1997

free call 1800 066 077

www.fisherbiotec.com

For all your life science research requirements

# Our latest SPECIAL OFFERS



[VIEW OUR ENTIRE PRODUCT RANGE](#)



**\$45**  
ex-GST  
per pack  
(10 racks x 96 tips)

**\*REDUCED\***



<< [VIEW FLYER](#)

**Sorenson 10µl Extra-long sharp tip MultiGuard™ Barrier Tips**

## FINAL WAREHOUSE CLEARANCE

*The most accurate sample delivery tip in the industry\**

- Ideal for volumes of less than 1µl
- Extended length, fine point
- Universal fit for all research-grade pipettors



Offer valid until stocks sold out

## Electroporation – Cuvettes Plus

*Gamma irradiated for guaranteed sterility!*

- 3 sizes available – 1 mm, 2 mm, and 4 mm
- Compatible with most popular brands of Electroporators
- Cuvettes can obtain high field strengths up to 25.0 kV/cm
- Disposable, simply discard carefully after use



[VIEW FLYER >>](#)



*Great value!*

Contact 1800 066 077 for pricing

## RNase and DNA Inactivation Spray – EliminatoR

*Convenient and ready-to-use!*

Use EliminatoR spray before performing PCR or working with RNase and DNA

Fast and easy decontamination to maintain an RNase- and DNA-free work area

Handy refill bottle available - \$66 (ex-GST)

[VIEW FLYER >>](#)



**\$70**  
ex-GST  
500ml bottle



## SERANA Foetal Bovine Serum

Australian origin and manufacture

- Certified free from all animal diseases
- Suitable for cell and tissue culture
- Triple 0.1µm
- Sterile filtered



Serum expiry date: May 2024

[VIEW FLYER >>](#)

**\$590**  
ex-GST  
500ml bottle



Offer valid while stocks last

Terms and Conditions

All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Fisher Biotec reserves the right to withdraw offers at any time.

## COVID-19 Virus Detection



## HighPrep™ Viral DNA/RNA Kits

**RAPID AND RELIABLE MAGNETIC BEAD BASED ISOLATION OF VIRAL NUCLEIC ACIDS**

- Efficient purification of viral nucleic acids
- Adaptable to various automated liquid handlers
- No toxic organic solvents



<< [FIND OUR MORE](#)

Contact **1800 066 077** for pricing


## iSWAB™ Microbiome Collection Kits

- Safer handling & fewer bio-hazard precautions
- Viral RNA stable at room temp. for up to 21 days
- Any swab can be used

Validated by the CDC for collection of Coronavirus samples, due to its property of inactivating the Virus and not requiring cold chain shipment



## COVID-19 Virus Detection



<< [FIND OUT MORE](#)

Contact **1800 066 077** for pricing

Contact **1800 066 077** for pricing



Available in 25 mL bottle as a 4X concentrate

'Unique Cell Culture Reagent'



<< [MORE INFORMATION](#)

## Happy Cell™ Advanced Suspension Medium (ASM)

*A 3D Culture and Cell Suspension Medium*

Happy Cell® ASM is:

- Very cost effective and easy-to-use
- Low viscosity and density
- Quick and easy sample recovery
- Formulated for automated liquid handling
- Use in any type of microplate assay

**NEW**

Release!



Contact **1800 066 077** for pricing

## Precision Plus Series Advanced BenchTop Pipettor

*'Accelerate your liquid transfers'*

**Available in two volume ranges: 25µl & 300µl**

- 96/384 semi-automated pipettor
- For use with any Universal tips
- Compact, easily fits into a fume hood
- 6 pipetting speeds

[VIEW THE FLYER >>](#)



## Axygen® Storage Boxes and Microcentrifuge Tube Racks

*Ideal for sample storage*

**Tube Racks –**

- Step design allows easy access to tubes
- 2 rows of 10 (20 x 1.5 to 2.0 mL total)



**Storage Boxes –**

- Freezing at -80°C and autoclaving up to 121°C
- Alpha-numeric sample identification
- 100 x 1.5 to 2.0 mL tube capacity

[FIND OUT MORE >>](#)

**NEW**  
Release  
from Axxygen®



**For more information go to [www.fisherbiotec.com](http://www.fisherbiotec.com)**

Terms and Conditions

All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Fisher Biotech reserves the right to withdraw offers at any time.