

# SPM, double-sided adhesive mat

## DESCRIPTION

Convenient alternative to traditional steel holders, an easy way to fix tubes, plates, flasks and other laboratory wares on platforms for cultivation, incubation and mixing.

SPM is compatible with PP-4 platform, that fits both on PSU-10i orbital shaker and in ES-20 Shaker-Incubator.

Made of polyurethane with adhesive, simple to clean and durable, able to withstand up to 1,000 times placement/removal or 12 months of use. Additional information about temperature, working volume and speed limitations available in user manual.



## SPECIFICATIONS

Colour	transparent
Duration of use	up to 1,000 times placement/removal or 12 months
Temperature range	+4C° ... +80C°
Speed range	0–300 rpm
Overall dimensions (W×D×H)	210 × 210 × 3 mm (double sided 1.5 mm PU with PET adhesive)

## CAT. NUMBER


BS-010111-BK SPM, Double-sided adhesive mat

## ACCESSORIES



**ES-20**  
Without platform  
Orbital Shaker-Incubator

The ES-20 is a compact bench-top Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program.



**PP-4**  
BS-010108-BK  
platform

Flat platform with non-slip silicone mat for Petri dishes, culture flasks, agglutination cards.



**PSU-10i**  
Without platform  
Orbital Shaker

Shaker PPSSUU-1100ii provides regulated orbital motion of the platform and is designed for use both in small specialized biotechnological laboratories and in large multidisciplinary laboratories.