

free call 1800 066 077

www.fisherbiotec.com

For all your life science requirements



# SERANA

## Foetal Bovine Serum



**REDUCED**  
from \$590

now **\$550**  
ex-GST  
(500 ml)



**EXPIRY DATE:**  
**MAY 2024**

Catalog No FBS-AU-015

### SPECIFICATIONS

Foetal Bovine Serum



pH	6.8–8.2
Osmolality	280–340 mOsm/kg
Endotoxin	< 10 EU/ml
Total Protein	3.0–4.5 g/dl
Albumin	1.4–3.4 g/dl
IgG	< 400 µg/ml
Free Haemoglobin	< 25 mg/dl
Mycoplasma	Not detected
Virus & Antibodies	BHV-1, PI-3, BVDV-1,2
Sterility	Tested
Long Term Storage	< -15°C

[View the Certificate of Analysis](#)

SERANA Serum is of AUSTRALIAN origin and manufacture

- Certified free from all animal diseases
- Triple 0.1µm
- Sterile Filtered
- Suitable for cell and tissue culture

SERANA'S Foetal Bovine Serum is sourced from only the highest quality raw material sources in Australia. SERANA have extensive industry experience manufacturing to the highest of standards. SERANA Foetal Bovine Serum is affordable and Quality Assured, with extensive Certificate of Analysis.

**BUY NOW**

**OFFER VALID WHILE STOCKS LAST**

Terms and Conditions

All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Fisher Biotec reserves the right to withdraw offers at any time.