



free call 1800 066 077

www.fisherbiotec.com

For all your life science requirements

# CelCulture<sup>®</sup> CO<sub>2</sub> Incubator




**ESCO**  
WORLD CLASS. WORLDWIDE.

**\$6,900**  
(ex-GST)



## CO<sub>2</sub> Incubator FOR CELL / TISSUE CULTURE

- Features a 170 litre full stainless steel chamber with easy-clean rounded corners.
- Dual decon systems.
- Automatic with internal ULPA filter and user activated 90°C moist heat.
- Latest generation dual wavelength, high heat, maintenance free IR sensor.
- Glass inner door and 30mm access port.
- On board data-logging and two year warranty.



Order 2 units and pay \$6,600 ea.  
(includes roller base stand)

### ORDERING INFORMATION

Catalog No	Description	Qty.	Special Price (ex-GST)
2170002	CelCulture <sup>®</sup> Incubator, 170L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz	1	\$6,900 ea.
2170002 (with roller-base)	CelCulture <sup>®</sup> Incubator, 170L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz, includes roller base stand*	2 (buy 2 units*)	\$6,600 ea.

Please quote Catalog No when ordering

**BUY NOW**

This offer expires on 30 June 2021

Terms and Conditions

All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Fisher Biotec reserves the right to withdraw offers at any time.