

FB Maestrogen DNA Ladders

fisherbiotec
AUSTRALIA


LOW PRICE!

Accurate size and mass estimations


Only \$85 (ex GST)

DNA Ladders comprising of fragments of specific base pair length are available for estimation of molecular weights from 100bp to 10Kb.

- For easy determination of band size, these DNA Ladders include bands that are two to three times more intense.



AccuRuler 100 bp Plus DNA Ladder



AccuRuler 1Kb Plus DNA Ladder

Catalog No.	Description	Quantity	Price (ex GST)
02002-500	AccuRuler 100 bp Plus DNA - 100 bp to 3,000 bp. The 500 bp and 1,500 bp bands have increased intensity	50µg / 500µl	\$85
02004-500	AccuRuler 1Kb Plus DNA Ladder - 100 bp to 10 kb. The 1 kb and 3kb bands have increased intensity	50µg / 500µl	\$85

