



free call 1800 066 077

www.fisherbiotec.com

For all your life science requirements

# CelCulture<sup>®</sup> CO<sub>2</sub> Incubator



**ESCO**  
WORLD CLASS. WORLDWIDE.

## FOR CELL / TISSUE CULTURE

**\$6,900**  
(ex-GST)



- ✓ 170L full stainless steel chamber with easy-clean rounded corners.
- ✓ Dual decontamination systems –
  - Internal ULPA filter to ISO Class 5 cleanliness
  - User activated automatic 90°C moist heat cycle
- ✓ Latest generation IR CO<sub>2</sub> Sensor –
  - Dual wavelength
  - Auto zeroing
  - Drift-free
  - High heat stable
  - Maintenance free
- ✓ Glass inner door and 30mm access port
- ✓ On board data-logging and 2 year warranty



Order 2 units and pay \$6,600 ea.  
(includes roller base stand)

### ORDERING INFORMATION

Catalog №	Description	Qty.	Special Price (ex-GST)
2170002	CelCulture <sup>®</sup> Incubator, 170L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz	1	\$6,900 ea.
2170002 (with roller-base)	CelCulture <sup>®</sup> Incubator, 170L, IR Sensor, CO <sub>2</sub> Control, ULPA, Moist Heat Decon, 230 VAC, 50/60 Hz, includes roller base stand*	2 (buy 2 units*)	\$6,600 ea.

Please quote Catalog № when ordering

**BUY NOW**

This offer expires on 30 Sep 2021

Terms and Conditions

All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Fisher Biotec reserves the right to withdraw offers at any time.