

BTX® Cuvettes Plus

free call **1800 066 077**

ELECTROPORATION CUVETTES PLUS

For use in electroporation and electrofusion of bacteria, yeast, insect, plant and mammalian cells

FEATURES:

- Cuvettes and safety stand protect both user and sample
- Compatible with most popular commercially available electroporators
- Cuvettes can obtain high field strengths up to 25.0kV/cm
- Gamma irradiated for guaranteed sterility
- Available in three electrode gap sizes, 1mm, 2mm and 4mm
- Disposable, simply discard carefully after use.



Ordering information

Item No	Description
45-0124	Cuvettes Plus, 1 mm gap, 90 µl, sterile, individually wrapped, includes Transferpette, Pkg/50, Grey
45-0125	Cuvettes Plus, 2 mm gap, 400 µl, sterile, individually Wrapped, includes Transferpette, Pkg/50, Blue
45-0126	Cuvettes Plus, 4 mm gap, 800 µl, sterile, individually Wrapped, includes Transferpette, Pkg/50, Yellow

ORDER NOW

Please contact us for pricing on 1800 066 077 or email info@fisherbiotec.com