



free call **1800 066 077**

www.fisherbiotec.com

For all your life science research requirements

Promotion

PROTEIN RESEARCH TOOLS

abcam




Lightning-Link® Antibody, Protein & Peptide Conjugation kits

Up to
30%
off RRP

OFFER EXPIRES 31 OCTOBER 2021

- Quickly conjugate the included label to your antibody, protein or peptide in as little as 20 minutes
- 100% antibody recovery in a wide range of applications, including FACS, ELISA, IHC, WB and Lateral flow
- Economical, can be used with just 10µg of antibody

Quote promo code: **50OFF-X5ZKB**



DyLight® 550 Conjugation Kit (Fast) - Lightning-Link® (ab201800)

Protein Ladders Prestained & Unstained

Great
Value!


OFFER EXPIRES 31 OCTOBER 2021

- Quick, ready-to-use, no heating or dilution
- Strong band intensity for high transfer efficiency
- Sharp bands for accurate MW determination

Description – Pre-stained Protein Ladders	Catalog Nº	Special Price*
Mid-range molecular weight (10 – 180 kDa)	ab116027	\$149.00
Broad molecular weight (10 – 245 kDa)	ab116028	\$158.00
Extra broad molecular weight (5 – 245 kDa)	ab116029	\$157.00
Extra broad molecular weight (6.5 – 270 kDa)	ab234592	\$154.00
Mid-range molecular weight (10 – 180 kDa)	ab234617	\$153.00
Unstained Protein Ladder (10 – 200 kDa)	ab234618	\$153.00

*Pricing is ex-GST

Quote promo code: **LADDER-X5ZKV**



SDS-PAGE - Prestained Protein Ladder – Broad molecular weight (10-245 kDa) (ab116028)

free call **1800 066 077** for pricing inquiries or email info@fisherbiotec.com

InstantBlue®

Ultra fast protein stain

NOW
20%
off RRP

OFFER EXPIRES 31 OCTOBER 2021

abcam

One step, ultra fast
InstantBlue®
Coomassie Protein Stain



- Ready-to-use Coomassie protein stain for polyacrylamide gels with results in 15 mins
- Remove SDS, fix and stain in ONE STEP
- No de-stain is needed
- As little as 5ng/band (BSA) of protein can be detected with low background
- Non-toxic – contains no methanol, suited for regular use and disposal

Description – Ultra Fast Protein Stain	Catalog №	Special Price*
InstantBlue Coomassie Protein Stain (1Litre)	ab11911	\$234.00

*Pricing is ex-GST




RunBlue®

Precast Gels

Offering superior rigidity & stability

- Sharper bands and batch-to-batch consistency
- Long shelf life, 10 times stronger than other brands
- Optimized buffers – running buffer offers low heat generation
- Compatible with all commonly used electrophoresis equipment



free call **1800 066 077** for pricing inquiries or email info@fisherbiotec.com

www.fisherbiotec.com