


# Combs for Vertical Systems

Cleaver Scientific offer the widest range of combs available of any gel tank manufacturer. They fit virtually every application from preparatory electrophoresis to high-throughput screening.

Available in four thicknesses and are color-coded:

- White - 1 mm supplied as standard
- Black - 0.7 mm for tightly resolved bands
- Red - 1.5 mm to maximize sample volume
- Blue - 2 mm to maximize sample volume.



omniPAGE Mini			
Catalog №	№ wells	Comb thickness	Well spacing
VS10-1-0.75	1	0.75 mm	Standard
VS10-5-0.75	5	0.75 mm	Standard
VS10-9-0.75	9	0.75 mm	Standard
VS10-10-0.75	10	0.75 mm	Standard
VS10-12-0.75	12	0.75 mm	Standard
VS10-16-0.75	16	0.75 mm	Standard
VS10-20-0.75	20	0.75 mm	Standard
VS10-8MC-0.75	8	0.75 mm	Multi-channel
VS10-1-1	1	1 mm	Standard
VS10-5-1	5	1 mm	Standard
VS10-9-1	9	1 mm	Standard
VS10-10-1	10	1 mm	Standard
VS10-12-1	12	1 mm	Standard
VS10-16-1	16	1 mm	Standard
VS10-20-1	20	1 mm	Standard
VS10-8MC-1	8	1 mm	Multichannel
VS10-1-1.5	1	1.5 mm	Standard
VS10-5-1.5	5	1.5 mm	Standard
VS10-9-1.5	9	1.5 mm	Standard
VS10-10-1.5	10	1.5 mm	Standard
VS10-12-1.5	12	1.5 mm	Standard
VS10-16-1.5	16	1.5 mm	Standard
VS10-20-1.5	20	1.5 mm	Standard
VS10-8MC-1.5	8	1.5 mm	Multichannel
VS10-1-2	1	2 mm	Standard
VS10-5-2	5	2 mm	Standard
VS10-9-2	9	2 mm	Standard
VS10-10-2	10	2 mm	Standard
VS10-12-2	12	2 mm	Standard
VS10-16-2	16	2 mm	Standard
VS10-20-2	20	2 mm	Standard
VS10-8MC-2	8	2 mm	Multichannel
omniPAGE Maxi Plus			
Catalog №	№ wells	Comb thickness	Well spacing
VS30-1-1	1	1 mm	Standard
VS30-2-1	2	1 mm	Standard
VS30-4-1	4	1 mm	Standard
VS30-75-1	75	1 mm	Standard
VS30-28MC-1	28	1 mm	Multichannel
VS30-56MC-1	56	1 mm	Multichannel

omniPAGE Mini Wide & omniPAGE Wave Maxi			
Catalog №	№ wells	Comb thickness	Well spacing
VS20-1-0.75	1	0.75 mm	Standard
VS20-5-0.75	5	0.75 mm	Standard
VS20-10-0.75	10	0.75 mm	Standard
VS20-20-0.75	20	0.75 mm	Standard
VS20-30-0.75	30	0.75 mm	Standard
VS20-36-0.75	36	0.75 mm	Standard
VS20-48-0.75	48	0.75 mm	Standard
VS20-18MC-0.75	18	0.75 mm	Multi-channel
VS20-1-1	1	1 mm	Standard
VS20-5-1	5	1 mm	Standard
VS20-10-1	10	1 mm	Standard
VS20-20-1	20	1 mm	Standard
VS20-30-1	30	1 mm	Standard
VS20-36-1	36	1 mm	Standard
VS20-48-1	48	1 mm	Standard
VS20-18MC-1	18	1 mm	Multichannel
VS20-1-1.5	1	1.5 mm	Standard
VS20-5-1.5	5	1.5 mm	Standard
VS20-10-1.5	10	1.5 mm	Standard
VS20-20-1.5	20	1.5 mm	Standard
VS20-30-1.5	30	1.5 mm	Standard
VS20-36-1.5	36	1.5 mm	Standard
VS20-48-1.5	48	1.5 mm	Standard
VS20-18MC-1.5	18	1.5 mm	Multichannel
VS20-1-2	1	2 mm	Standard
VS20-5-2	5	2 mm	Standard
VS20-10-2	10	2 mm	Standard
VS20-20-2	20	2 mm	Standard
VS20-30-2	30	2 mm	Standard
VS20-36-2	36	2 mm	Standard
VS20-48-2	48	2 mm	Standard
VS20-18MC-2	18	2 mm	Multichannel
omniPAGE Maxi Plus (contin.)			
Catalog №	№ wells	Comb thickness	Well spacing
VS30-1-1.5	1	1.5 mm	Standard
VS30-2-1.5	2	1.5 mm	Standard
VS30-4-1.5	4	1.5 mm	Standard
VS30-75-1.5	75	1.5 mm	Standard
VS30-28MC-1.5	28	1.5 mm	Multichannel
VS30-56MC-1.5	56	1.5 mm	Multichannel

Please contact us for pricing inquiries on 1800 066 077