



Fisher Biotec RNase and DNA EliminatoR Spray



www.fisherbiotec.com

For all your life science requirements

Fisher Biotec offers a comprehensive range of quality products at competitive prices. We specialise in the manufacture of leading-edge products for Molecular Biology, Genomics and Proteomics research.

Rnase and DNA Inactivation Spray - EliminatoR

Fisher Biotec DNA/RNA *EliminatoR* Inactivation Spray effectively eliminates DNase and DNA from laboratory glassware and plastic ware products without the need for autoclaving and eliminates time spent baking glassware.

EliminatoR should be used before performing PCR or working with DNA for a RNase and DNA-free work area. *EliminatoR* also eliminates the need for carcinogenic DEPC treatments.

The ready-to-use formula can be applied to surfaces to decontaminate, and wipe dry or rinse clean.

Catalog No	Description	Pack Size
RNX-500	Rnase/DNA EliminatoR, spray bottle	500 ml
RNX-500-R	Rnase/DNA EliminatoR, refill bottle	500 ml
RNX-1000-R	Rnase/DNA EliminatoR, refill bottle	1000 ml



ORDER NOW

For more information about our products – free call 1800 066 077