



BN-13/75
BS-010208-PK
adapter set

For rotors R-12/10, Ri-12/10
optional adapter set (12 pcs) -
BN-13/75 for vacutainers 2-5 ml



BN-13/100
BS-010208-QK
adapter set

For rotors R-12/10, Ri-12/10
optional adapter set (12 pcs) -
BN-13/100 for vacutainers 4-
8 ml



BN-16/100
BS-010208-RK
adapter set

For rotors R-12/10, Ri-12/10
optional adapter set (12 pcs) -
BN-16/100 for vacutainers 8-10
ml



Ri-24/10
BS-010221-LK
Rotor

Swing-out rotor for round
bottom centrifuge tube.



BN-13/75-24
BS-010221-RK
adapter set

For rotors R-24/10, Ri-24/10
optional adapter set (24 pcs) -
BN-13/75-24 for vacutainers 2-5
ml



BN-13/100-24
BS-010221-SK
adapter set

For rotors R-24/10, Ri-24/10
optional adapter set (24 pcs) -
BN-13/100-24 for vacutainers 4-
8 ml



BN-16/100-24
BS-010221-TK
adapter set

For rotors R-24/10, Ri-24/10
optional adapter set (24 pcs) -
BN-16/100-24 for vacutainers 8-
10 ml



Ri-12/15
BS-010221-KK
Rotor

Swing-out rotor for conical
bottom centrifuge tube.



Ri-2
BS-010221-GK
Rotor

Swing-out rotor for standard 96-
well microtitre plates, skirted
PCR plates, deepwell plates up
to 45mm.



AP-96
BS-010219-DK
adapter

2 adapters for 96-well semi-
skirted and non-skirted PCR
plates - made of Ertacetal® C.
Autoclavable



AP-384
BS-010219-EK
adapter

2 adapters for 384-well plates -
made of Ertacetal® C.
Autoclavable



Ri-24GC
BS-010221-NK
Rotor

Rotor for 24 Gel Cards for blood
group serology testing (Forward
Group, Reverse Group, RhD
Type and 3 cell antibody
screen)



RR-U
BS-010208-UK
rack for rotors

Rack for rotors