

## Labnet *Force Mini* Personal Centrifuge

The Labnet Force Mini Microcentrifuge is a personal centrifuge that is ideal for quick spin downs of microtubes and PCR Tubes. Rotors are quickly and easily interchanged. The centrifuge was designed to be extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpms. Press the lid release button to open the lid and the rotor quickly decelerates to stop for removing samples.

- Ideal for microfiltration and quick spins
- Complete with 8 x 2.0mL microtube rotor and 4 x 8 PCR strip rotor
- Quickly reaches maximum speed of 6,000 rpm/2,000 x g
- Transparent lid for easy sample viewing.



Force Mini Microcentrifuge

### SPECIFICATIONS

#### *Force Mini* Centrifuge

Maximum speed/RCF

6,000 rpm/2,000 x g

Tube capacity

8 x 1.5/2.0mL tubes; 4-8 x 0.2mL tube strips or 32 x 0.2mL tubes

| CAT NO.  | DESCRIPTION  |
|----------|--|
| SBC140-3 | Labnet Force Mini Microcentrifuge, with 1.5mL tube and strip tube rotors, 8 adapters for 0.4 and 0.5mL tubes, 230V |
| SBC140TR | Additional 8 x 1.5mL tube rotor for Force Mini   |
| C1205    | Adapters for 0.5/0.6mL tubes, pack of 6  |
| C1206    | Adapters for 0.25/0.4mL tubes, pack of 6   |
| C1222    | Adapters for 0.2mL tubes, pack of 6  |