



SPECIAL OFFER UV Air Purifiers

What are the benefits of UVR-M & UVR-Mi Air Purifiers?

The Biosan UV Air Purifiers are ideal for sanitising air in most indoor environments. They are safe to use and effective against common airborne diseases by sanitising the air and efficiently destroying disease-causing agents (e.g. Viruses, Micro-organisms) by UV radiation. ***Published research indicates that UV air disinfection may be an effective tool for preventing important respiratory viral diseases such as SARS-Coronavirus.***

- Full user protection from direct UV radiation
- Very low noise levels
- Simple to operate and maintain
- Digital control unit can track UV lamp operating time (UVR-Mi only)
- Easy to install – fix to a wall, or mount on the movable stand.



\$1,170.00

UVR-Mi

\$675.00

UVR-M

SPECIFICATIONS	UVR-M	UVR-Mi
UV radiation source: 25W Bactericidal (Philips)	1 lamp	2 lamps
UV radiation level	18 mW / cm ² / sec	36 mW / cm ² / sec
Air-flow productivity	21m ³ / hour	29m ³ / hour
Full user protection from direct UV light		Yes
Display	—	LCD
UV lamp operation indicator		Yes
UV lamp lifetime counter	No	Yes
Timer	—	1 min-24 hrs / non-stop
Automatic switch ON/OFF	—	Yes
Lamp fault detection	—	Yes

ORDERING

CATALOG №	DESCRIPTION	QTY.	PRICE (ex-GST)*
BS-040110-AA3	Biosan UVR-Mi Air Purifier / UV Cleaner-Recirculator, (2 lamps), 110 x 135 x 660 mm	1	\$1,170.00
BS-040105-AA3	Biosan UVR-M Air Purifier / UV Cleaner-Recirculator, (1 lamp), 110 x 135 x 660 mm	1	\$675.00
BS-040105-AK	Biosan UVR-S tripod stand, 510 x 510 x 1100 mm	1	\$189.00

ORDER NOW

Email info@fisherbiotec.com for more details or visit our [website](http://www.fisherbiotec.com)